# **Embedding yeast colonies for light and electron microscopy**

Sarah Piccirillo, Melissa G. White, Jeffrey C. Murphy, Douglas J. Law, and Saul M. Honigberg

# **Abstract**

The following is a modification of a Scherz R., 2001 method for colony embedding. The protocol described here is from:

Sarah Piccirillo, et. al. The Rim101p/PacC Pathway and Alkaline pH Regulate Pattern Formation in Yeast Colonies (2010)

Genetics 184:707-716; doi:10.1534/genetics.109.113480

Please see the full manuscript for additional details.

**Citation:** Sarah Piccirillo, Melissa G. White, Jeffrey C. Murphy, Douglas J. Law, and Saul M. Honigberg Embedding yeast colonies for light and electron microscopy. **protocols.io** 

dx.doi.org/10.17504/protocols.io.dw97h5

Published: 04 Dec 2015

#### **Guidelines**

To embed colonies for sectioning, we modified a previously described method (<u>Scherz et al. 2001</u>), with a major change being the substitution of Spurr's reagent as embedding medium.

#### **Original method:**

SCHERZ, R., V. SHINDER and D. ENGELBERG, 2001 <u>Anatomical analysis of Saccharomyces cerevisiae stalk-like structures reveals spatial organization and cell specialization</u>. J. Bacteriol. 183: 5402–5413.

#### **Protocol**

#### Step 1.

Incubate approximately 300 colonies on agar medium for the indicated time.

#### Step 2.

Remove an isolated colony (1-2 mm in diameter) and a small amount of the underlying agar medium, and place on a microscope slide

#### Step 3.

Place several drops of 2% agar (42°C) on a microscope slide, and immediately place the colony on the agar and then place several drops of agar on top of the colony and allow to solidify.

#### Step 4.

Trim the resulting agar block with a razor blade, and place in a 3.5 ml borosilicate screw-cap vial (Fisher 03-339-21B) containing 1.5ml 2%paraformaldehyde/2%glutaraldehyde.

#### NOTES

Saul Honigberg 30 Sep 2015

All subsequent incubations and washes use 1.5 -2.0 ml and are performed in the same vial.

#### Step 5.

Fix colonies by incubating for 7 days at 4°C

#### Step 6.

Wash #1: Wash agar blocks on ice by incubating for 15 minutes with 1.5 ml of 0.15M sodium cacodylate (pH 7.2)

© DURATION

00:15:00

#### Step 7.

Wash #2: Wash agar blocks on ice by incubating for 15 minutes with 1.5 ml of 0.15M sodium cacodylate (pH 7.2)

**O DURATION** 

00:15:00

#### Step 8.

Wash #3: Wash agar blocks on ice by incubating for 5 minutes with 1.5ml OS buffer (100 mM KH<sub>2</sub>PO<sub>4</sub>, 10 mM MqCl<sub>2</sub>, pH 6.0).

**O DURATION** 

00:05:00

#### Step 9.

Wash #4: Wash agar blocks on ice by incubating for 5 minutes with 1.5ml OS buffer (100 mM KH<sub>2</sub>PO<sub>4</sub>, 10 mM MgCl<sub>2</sub>, pH 6.0).

O DURATION

00:05:00

# Step 10.

If sections will be used for electron microscopy, add 1% OsO<sub>4</sub> in OS to vials to cover the agar blocks and incubate on ice in a chemical fume hood for 1 hr. Otherwise skip to step 13.

© DURATION

01:00:00

#### **Step 11.**

Dispose of 1% OsO<sub>4</sub> in hazardous waste.

WashA: Wash with 1.5ml OS buffer by incubating on ice for 10 minutes.

#### Step 12.

WashB: Wash with 1.5ml OS buffer by incubating on ice for 10 minutes.

## **Step 13.**

Add 1.5 ml OS and incubate overnight at 4°C.

© DURATION

18:00:00

# Step 14.

WashA: Wash blocks with 1.5ml cold water by incubating on ice for 10 minutes.

#### **Step 15.**

WashB: Wash blocks with 1.5ml cold water by incubating on ice for 10 minutes.

## **Step 16.**

Wash #1: Add 1.5ml cold 25% ethanol and incubate on ice for 10 minutes.

**O DURATION** 

00:10:00

## **Step 17.**

Wash #2: Add 1.5ml cold 50% ethanol and incubate on ice for 10 minutes.

**O DURATION** 

00:10:00

# Step 18.

Wash #3: Add 1.5ml cold 75% ethanol and incubate on ice for 10 minutes.

**O** DURATION

00:10:00

# Step 19.

Wash #4: Add 1.5ml cold 95% ethanol and incubate on ice for 10 minutes.

**O DURATION** 

00:10:00

## Step 20.

Wash #5: Add 1.5ml cold 100% ethanol and incubate on ice for 10 minutes.

© DURATION

00:10:00

#### Step 21.

Wash #6: Add 1.5ml cold 100% ethanol and incubate on ice for 10 minutes.

**O** DURATION

00:10:00

#### Step 22.

Remove ethanol and resuspend in 1.5ml cold 100% ethanol. Leave the blocks overnight at 4°C in 100% ethanol.

**O DURATION** 

18:00:00

#### Spurr's treatment

## Step 23.

Make Spurr's reagent by stirring slowly under chemical fume hood 5 grams ERL4221, 4 grams DER736 and 13 grams NSA for 20 minutes(Electron Microscopy Sciences).

#### Spurr's treatment

# Step 24.

Add 0.15 grams DMAE(Electron Microscopy Sciences) and stir for 20 minutes. De-gas for 1-2 hrs.

**O DURATION** 

00:15:00

# Spurr's treatment

#### **Step 25.**

Wash blocks with 1.5ml unopened room temp 100% ethanol by incubating at room temp for 10 minutes. Repeat 4 more times.

**O DURATION** 

00:30:00

#### Spurr's treatment

## Step 26.

Treatment#1: Remove ethanol and add 1.5ml of a 2:1 ratio of 100% ethanol:Spurr's reagent.

© DURATION

00:15:00

# Spurr's treatment

#### **Step 27.**

Treatment#1: Rotate vial for 15 minutes on wheel at room temperature.

**O DURATION** 

00:30:00

## Spurr's treatment

## Step 28.

Treatment#1: Allow to stand for 30 minutes at room temperature.

© DURATION

00:15:00

#### Spurr's treatment

## Step 29.

Treatment#2: Remove 100%ethanol:Spurr's reagent and add 1.5ml of a 2:1 ratio of 100% ethanol:Spurr's reagent. Rotate vial for 15 minutes on wheel at room temperature.

**O** DURATION

00:30:00

## Spurr's treatment

# Step 30.

Treatment#2: Allow to stand for 30 minutes at room temperature.

© DURATION

00:15:00

#### Spurr's treatment

#### **Step 31.**

Treatment #3: Remove 100% ethanol: Spurr's reagent and add 1.5ml of a 2:1 ratio of 100% ethanol: Spurr's reagent. Rotate vial for 15 minutes on wheel at room temperature.

**O DURATION** 

00:30:00

## Spurr's treatment

#### **Step 32.**

Treatment #3: Allow to stand for 30 minutes at room temperature.

© DURATION

00:15:00

#### Spurr's treatment

#### Step 33.

Treatment #4: Remove 100% ethanol: Spurr's reagent and add 1.5ml of a 1:1 ratio of 100%ethanol: Spurr's reagent. Rotate vials for 15 minutes on wheel at room temperature.

**O DURATION** 

00:30:00

# Step 34.

Treatment #4: Allow to stand for 30 minutes at room temperature.

© DURATION

04:00:00

#### Step 35.

Treatment#5: Remove 100% ethanol: Spurr's reagent and add 1.5ml of a 1:1 ratio of 100% ethanol: Spurr's reagent. Rotate vials for 15 minutes on wheel at room temperature.

© DURATION

18:00:00

#### Step 36.

Treatment#5: Allow to stand for 30 minutes at room temperature.

**O DURATION** 

18:00:00

#### Step 37.

Remove 100% ethanol:Spurr's reagent and add 1.5ml Spurr's reagent to vial. Incubate for 4hrs at room temperature.

© DURATION

04:00:00

# Step 38.

Replace with 1.5ml Spurr's reagent and rotate the vial overnight on the wheel at room temperature.

© DURATION

18:00:00

#### Step 39.

Replace with 1.5ml Spurr's reagent and rotate the vial until late afternoon on the wheel at room temperature.

# Step 40.

Remove the Spurr's reagent and replace with 1.5ml freshly made Spurr's reagent. Rotate the vial overnight on the wheel at room temperature.

**O DURATION** 

00:00:15

# **Step 41.**

The next day replace with 1.5ml Spurr's reagent and rotate until the following day on the wheel at room temperature.

#### **Step 42.**

The next day replace with 1.5ml Spurr's reagent and rotate until the following day on the wheel at room temperature.

#### Step 43.

Place each agar block in a mold with 0.2ml Spurr's reagent and incubate at 60°C for four hours.

#### **Step 44.**

Top off the molds with Spurr's reagent and incubate at 60°C for 3 days.

#### Step 45.

Collect sections(0.5u) from the central region of the colony in a drop of dH20 on a glass slide. Dry slide on a 52°C heat block.

#### Step 46.

Stain sections with 1% toluidine blue,1%Sodium Borate for 5-15 seconds.

#### Step 47.

Wash slide under a stream of dH20. Dry on heat block. Cover in Permount(Fisher SP15-100). Examine by light microscopy.