



Oct 16, 2019

DNase I treatment in solution (after RNA extraction) [↗](#)

[Ben Kuipers](#)¹¹Wageningen University[1](#) *Works for me* dx.doi.org/10.17504/protocols.io.75ghq3w[iGEM Wageningen 2019](#) [Ben Kuipers](#) 

ABSTRACT

This is a protocol used for the Dnase treatment after RNA isolation with the RNeasy mini kit. Protocol is based on the RNase-Free DNase Set and the RNeasy Mini Handbook 10/2013, Qiagen 'Appendix E: DNase Digestion of RNA before RNA Cleanup'.

EXTERNAL LINK

<https://www.qiagen.com/ie/products/discovery-and-translational-research/dna-rna-purification/rna-purification/total-rna/rneasy-mini-kit/#resources>

- 1 Check the concentrations of DNA and RNA with Qubit.
- 2 DNase I Stock solution was made from the Promega Maxwell 16 LEV SimplyRNA kit DNase by adding 275 µL nuclease free water to the lyophilized DNase I.
- 3 Prepare a mastermix of 5 µL Buffer RDD (RNase-Free DNase Set) + 1.25 µL DNase I stock (Promega Maxwell 16 LEV simplyRNA kit) + eluted sample. Make the volume up to 50 µL with RNase-free water (RNase-Free DNase Set).
- 4 Incubate the mixture on room temperature for 60 minutes.
- 5 Check the concentrations of DNA and RNA with Qubit again.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited