

# Simple & rapid genotyping of marine microeukaryotes

Jernej Turnsek

## Abstract

This is an effective protocol for genotyping diatom exconjugants which we believe is readily applicable to any other genetically manipulated marine microeukaryote.

**Citation:** Jernej Turnsek Simple & rapid genotyping of marine microeukaryotes. **protocols.io**

dx.doi.org/10.17504/protocols.io.jcdcis6

**Published:** 12 Aug 2017

## Guidelines

1. We find [Phire Plant Direct PCR Master Mix](#) to be a top choice for direct genotyping of diatom exconjugants.
2. Review protocol details [here](#).
3. Note that we regularly run 10 µL reactions as opposed to 20 µL which is the minimum volume recommended by the manufacturer.

## Before start

1. Your exconjugant colonies should be growing as small - ~300 µL - liquid cultures in the presence of appropriate antibiotics. A 48-well plate works great for this purpose.
2. Adjust primer T<sub>m</sub> using [this T<sub>m</sub> Calculator](#).

## Protocol

Prepare PCR master mix.

### Step 1.

--> volumes are per reaction

--> x (number of strains + 2)

☐ AMOUNT

4 µl Additional info: MQ

☐ AMOUNT

0.5 µl Additional info: 10 µM forward primer

☐ AMOUNT

0.5 µl Additional info: 10 µM reverse primer

☐ AMOUNT

5 µl Additional info: Phire Plant Direct PCR Master Mix

Aliquot into PCR tubes and add 0.5 µL of strains you wish to genotype.

### Step 2.

--> Include appropriate positive (if available) and negative controls.

--> Include WT or knockout strain you used for transformation as a template in one of the reactions.

Passage genotype+ strains.

### Step 3.

--> 1:100 or 1:50 dilution into 10 mL fresh medium is a good place to start.

Proceed with phenotypic characterization.

### Step 4.