

# Immunohistochemical staining

Xiao Bing Tang,Li Huan,Jin Zhang,Wei Lin Wang,Zheng Wei Yuan,Yu Zuo Bai

## Abstract

**Citation:** Xiao Bing Tang,Li Huan,Jin Zhang,Wei Lin Wang,Zheng Wei Yuan,Yu Zuo Bai Immunohistochemical staining. **protocols.io**

[dx.doi.org/10.17504/protocols.io.kujcwun](https://dx.doi.org/10.17504/protocols.io.kujcwun)

**Published:** 17 Nov 2017

## Protocol

### Step 1.

Slides are deparaffinized by immersion in xylene 2 times for 2 hr.

### Step 2.

Slides are rehydrated by immersion in 100% ethanol (10 min.), in 95% ethanol (10 min), in 85% ethanol (10min). Slides are washed with PBS (3 times, 3 min).

### Step 3.

Antigen retrieval is performed by heating the slides in 10 mM citrate buffer (pH 6.0) at 98°C for 10 min in a microwave oven. Slides are washed with PBS (3 times, 3 min).

### Step 4.

The endogenous peroxidase activity was blocked by incubation in 3% H<sub>2</sub>O<sub>2</sub> for 20 min. Slides are washed with PBS (3 times, 3 min).

### Step 5.

Antibody incubations were performed in phosphate-buffered saline (PBS) supplemented with 10% goat serum for 20 min at room temperature.

### Step 6.

Slides are incubated 16hr at 4°C in a humidifying box with primary antibody. Negative controls are made with PBS alone.

### Step 7.

Slides are washed with PBS (3 times, 3 min).

### Step 8.

Slides are incubated with biotinylated secondary antibody for 30 min in a humidifying box at room temperature.

#### **Step 9.**

Slides are washed with PBS (3 times, 3 min).

#### **Step 10.**

Slides are incubated with horseradish peroxidase (HRP)-conjugated secondary antibody (Santa Cruz Bio-technology) for 20 min in a humidifying box at room temperature.

#### **Step 11.**

Slides are washed with PBS (3 times, 3 min).

#### **Step 12.**

Signals were visualized by using 3 '3Pdiaminobenzidine (DAB; Sigma, UK) for 1 minute and terminated by incubated in distilled water.

#### **Step 13.**

Slides are washed with distilled water and counterstained with hematoxylin for 10 seconds.

#### **Step 14.**

Slides are washed with running water for 6 min.

#### **Step 15.**

Slides are dehydrated by immersion in 70% ethanol (1 min), in 95% ethanol (5 min), in 100% ethanol (2 times, 5 min) and in xylene (2 times, 20 min).

#### **Step 16.**

Cover slides are mounted with gum.