



## Prepare 2L of 2X freezer buffer to freeze down *C. elegans*

Cancer Research UK / Wellcome Gurdon Institute media kitchen<sup>1</sup>

<sup>1</sup>Cancer Research UK / Wellcome Gurdon Institute

[dx.doi.org/10.17504/protocols.io.3nhgmb6](https://doi.org/10.17504/protocols.io.3nhgmb6)

Working



Cristian Riccio ⚡

### ABSTRACT

This protocol accompanies you in the preparation of 2 litres of 2X freezer buffer, which is the buffer used to freeze down *C. elegans* strains for long-term storage

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
double distilled water (ddH2O)		
Sodium chloride meets analytical specification of Ph.Eur Fisher Chemical	S/3160/65	Fisher Scientific
Magnesium Sulfate Heptahydrate Certified AR for Analysis Fisher Chemical	M/1050/53	Fisher Scientific
Potassium Dihydrogen Orthophosphate Certified AR for Analysis Fisher Chemical	P/4800/53	Fisher Scientific
SYCHEM autoclave	View	Syschem
Glycerol 99 % Certified AR for Analysis Fisher Chemical	G/0650/17	Fisher Scientific

### 1 Freezer\_buffer.xls

### 2

Ingredients	Quantity	
NaCl 1M		200ml
1M KH <sub>2</sub> PO <sub>4</sub>	pH 6	100ml
Glycerol		600ml
Double distilled H <sub>2</sub> O	up to 2L	
		Autoclave
0.1M MgSO <sub>4</sub>	300 microlitres	
		per 100ml bottle

### 3

1	Combine NaCl, KH <sub>2</sub> PO <sub>4</sub> and Glycerol solutions with approx 900ml double distilled H <sub>2</sub> O in 2L measuring cylinder on magnetic stirrer.
2	Make up to volume using double distilled H <sub>2</sub> O
3	Aliquot into 20x100ml media bottles, label, date and autoclave.
4	After autoclaving and the solution is cold, add 300 microlitres of 0.1m MgSO <sub>4</sub> to each 100ml bottle using a Gilson in the flow hood.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited