

Sep 27, 2019

Priming and loading a MinION flowcell

[Josh Quick](#)¹¹University of Birmingham

1

Works for me

[dx.doi.org/10.17504/protocols.io.7q5hmy6](https://doi.org/10.17504/protocols.io.7q5hmy6)

ARTIC



Josh Quick ⚡ 🌞 🌱

- 1 Thaw the following reagents at room temperature before placing on ice:

Sequencing buffer (SQB)

Loading beads (LB)


Flush buffer (FLB)

Flush tether (FLT)

- 2 Add  30 µl FLT to the FLB tube and mix well by vortexing.


- 3 If required place a new MinION flowcell onto the MinION by flipping open the lip and pushing one end of the flowcell under the clip and pushing down gently.

- 4 Rotate the inlet port cover clockwise by 90° so that the priming port is visible.

- 5 Take a P1000 pipette and tip and set the volume to  800 µl . Place the tip in the inlet port and holding perpendicularly to the plane of the flowcell remove any air from the inlet port by turning the volume dial anti-clockwise.








Be careful not to remove so much volume that air is introduced onto the rectangular array via the outlet.

- 6 Load  800 µl of FLB (plus FLT) into the flow cell via the inlet port, dispense slowly and smoothly trying to avoid the introduction of any air bubbles.

- 7 Wait for  00:05:00 .

- 8 Gently lift the SpotON cover to open the SpotON port.


- 9 Load another  **200 µl** of FLB (plus FLT) into the flow cell via the inlet port, this will initiate a siphon at the SpotON port to allow you to load the library dilution.
- 10 In a new tube prepare the library dilution for sequencing:

Component	Volume
SQB	 37.5 µl
LB	 25.5 µl
Final library	 12 µl
Total	 75 µl



Mix LB immediately before use as they settle quickly.

Dilute library in EB if required.

- 11 Mix the prepared library gently by pipetting up and down just prior to loading.
- 12 Add the  **75 µl** library dilution to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop siphons into the port before adding the next.
- 13 Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the inlet port and close the MinION lid.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited