Flagellar transformation of Naegleria gruberi

Anastasios Tsaousis

Abstract

Citation: Anastasios Tsaousis Flagellar transformation of Naegleria gruberi. protocols.io

dx.doi.org/10.17504/protocols.io.hnhb5b6

Published: 31 Jul 2017

Protocol

Remove media (supernatant) from the flask

Step 1.

Wash with 2 mM Tris, pH7.6

Step 2.

Remove supernatant from the flask and add again 2 mM Tris, pH7.6

Step 3.

Incubate without shaking for at least 2 hours at 30 oC

Step 4.

Check for flagellation under a light microscope microscope

Step 5.

Use x10 or x20 objectives