



Jul 08, 2019

03 Ligation

TJUSLS China¹¹Tianjin University

Working

[dx.doi.org/10.17504/protocols.io.49igz4e](https://doi.org/10.17504/protocols.io.49igz4e)TJUSLS China
Tianjin University

MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

ddH₂O

Water bath or heat block at 37°C

5×T4 DNA Ligase Buffer

/

T4 DNA Ligase

/

- 1 Mix together:
5×T4 DNA Ligase Buffer 2μl
T4 DNA Ligase 0.5–1μl
Vector as required
Insert as required
ddH₂O Add to 10μl

- 2 Heat in water bath or heat block at 22°C for 30min.

🔥 22 °C

🕒 00:30:00



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited