



Jun 11,
2019

Poly-d lysine coating slides

Grace Burgin¹

¹Regev Lab

Working

[dx.doi.org/10.17504/protocols.io.32cgqaw](https://doi.org/10.17504/protocols.io.32cgqaw)



Grace Burgin



ABSTRACT

coating standard blank slides with poly-d lysine to increase adhesion of tissue

- 1 Make .3mg/ml in 0.1M borate buffer [→ add 16.67 ml 0.1M borate to 5mg bottle of poly-D lysine lyophilized powder]
- 2 coat the slides once by dipping
- 3 Dry
- 4 Rinse twice with dH₂O
- 5 Dry



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited