

Ki-67 Staining Protocol

Kelsey Miller

Abstract

Citation: Kelsey Miller Ki-67 Staining Protocol. [protocols.io](https://doi.org/10.17504/protocols.io.fabbian)

[dx.doi.org/10.17504/protocols.io.fabbian](https://doi.org/10.17504/protocols.io.fabbian)

Published: 08 Jul 2016

Protocol

Step 1.

Prepare 70% ethanol and chill at -20°C.

Step 2.

Prepare target cells of interest and wash 2X with PBS by centrifuge at 350xg for 5 minutes.

 [DURATION](#)

00:05:00

Step 3.

Discard supernatant and loosen the cell pellet by vortexing.

Step 4.

Add 3 ml cold 70% ethanol drop by drop to the cell pellet while vortexing.

Step 5.

Continue vortexing for 30 seconds and then incubate at -20°C for 1 hour.

 [DURATION](#)

01:00:00

Step 6.

Wash 3X with BioLegend's Cell Staining Buffer (Cat. No. 420201) and then resuspend the cells at the concentration of $0.5-10 \times 10^6/\text{ml}$.

Step 7.

Mix 100 μl cell suspension with proper fluorochrome-conjugated Ki-67 antibody and incubate at room temperature in the dark for 30 minutes.

 [DURATION](#)

00:30:00

Step 8.

Wash 2X with BioLegend's Cell Staining Buffer (Cat. No. 420201) and then resuspend in 0.5 ml cell staining buffer for flow cytometric analysis.