



Oct 15, 2019

## Antibiotic Stocks

Alba Balletbó<sup>1</sup><sup>1</sup>Wageningen University**1** Works for me dx.doi.org/10.17504/protocols.io.7krhkv6

iGEM Wageningen 2019



Alba Balletbó ⚡

### MATERIALS TEXT

Ampicillin  
Chloramphenicol  
Kanamycin  
Spectinomycin  
Tetracycline

- 1 Unless otherwise indicated, the antibiotic powder can be dissolved in dH<sub>2</sub>O. Addgene recommends making 1000X stock solutions and storing aliquots at -20°C.

Antibiotic	Recommended Stock Concentration	Recommended Working Concentration
Ampicillin	100 mg/mL	100 µg/mL
Chloramphenicol	25 mg/mL(dissolve in EtOH)	25 µg/mL
Gentamycin	10 mg/mL	10 mg/mL
Kanamycin	50 mg/mL	50 µg/mL
Spectinomycin/Streptomycin	50 mg/mL	50 µg/mL
Tetracycline*	10 mg/mL	10 µg/mL

\* Light sensitive, Mg<sup>2+</sup> - inhibitor - do not use with minimal media.

**'Note':** All are dissolved in MQ except Chloramphenicol, which is dissolved in ethanol. All antibiotic mixes were filter-sterilized before use.

- 2 Weigh the antibiotic powders needed for 10 mL of Recommended Stock Concentration Solution.
- 3 Add the weighed powder to the MQ or ethanol solution. Mix/vortex vigorously until the powder is completely dissolved.
- 4 Filter sterilize and aliquot into 1.5 mL tubes.
- 5 To use, dilute your antibiotic into your LB medium at 1:1,000. For example, to make 100 mL of LB/ampicillin growth media, add 100 µL of a 100 mg/mL ampicillin stock (1000X stock) to 100 mL of LB.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited