



## Golden Gate Cloning LVL 0 Version 2

Daniel Marchal<sup>1</sup>

<sup>1</sup>iGEM Team Marburg 2018

[dx.doi.org/10.17504/protocols.io.uvhew36](https://doi.org/10.17504/protocols.io.uvhew36)

Version 2

Oct 21, 2018

iGEM Marburg 2018

Working

Daniel Marchal

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS

NAME	CATALOG #	VENDOR
T7 DNA Ligase - 100,000 units	M0318S	New England Biolabs
nuclease free water		Contributed by users
Esp3I	R0734S	New England Biolabs
10X NEB T4 DNA ligase buffer		New England Biolabs


### Reaction Setup on ice


1. Add 1  $\mu$ L of 70 ng Template DNA. 0 °C
2. Add three fold excess of PCR-Fragment. 0 °C
3. Add 0.5  $\mu$ L Esp3I. 0 °C
4. Add 0.5  $\mu$ L T7-Ligase. 0 °C
5. Add 1  $\mu$ L T4-Ligase Buffer. 0 °C
6. Fill with Nuclease-free water to 10  $\mu$ L. 0 °C


## 7 Thermocycling conditions:

### Thermocycling conditions

8 30 Cycles of 2min 98°C / 5min 16°C

9 30 min. 37°C.  37 °C

10 10 min. 80°C.  80 °C

11 Hold 20°C.  20 °C



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited