

Protocol for drug sensitivity assay with *S. robusta*

Dan Needleman

Abstract

Citation: Dan Needleman Protocol for drug sensitivity assay with *S. robusta*. **protocols.io**

dx.doi.org/10.17504/protocols.io.g9ebz3e

Published: 23 Feb 2017

Protocol

Developed by Aaron Turkewitz and Lev Tsypin

Step 1.

- 1) Measure density of floating cells in culture (D6 strain):
 - a. Grow 100 mL volume of cells in a 175 cm² flask for three to four days, until there is a substantial number of healthy, floating cells
 - b. Stand flask vertically and take a volume of the floating culture (so as to not disturb the adherent cells) to measure cell density
 - c. Count floating cells (you may need to concentrate them, prior to counting, e.g., in a haemocytometer)
- 2) Distribute cells on a petri dish:
 - a. The target concentration for distributing cells is 400-500 floating cells/mL, so dilute as necessary. The volume needed for a plate is 2 mL, including the volume of any added drugs.
 - b. Transfer 2 mL, containing the desired number of cells, to an Eppendorf tube.
 - i. Centrifuge using a personal mini-centrifuge for two minutes. At this speed (which on our model is 6000 rpm, 1400 x g), many cells remain in the supernatant. Those are the cells to distribute in drops. Do not disturb cells in the pellet.

Developed by Aaron Turkewitz and Lev Tsypin

Step 2.