



Oct 16, 2019

## Chitin binding assay using chitin magnetic beads (NEB)

Robert Hoofman<sup>1</sup><sup>1</sup>Wageningen University

In Development

dx.doi.org/10.17504/protocols.io.8anhdsde

iGEM Wageningen 2019

Robert Hoofman  
Wageningen University

## MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Chitin Magnetic Beads - 5 ml

E8036S

New England Biolabs

## STEPS MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Chitin Magnetic Beads - 5 ml

E8036S

New England Biolabs

## 1 Make the following buffer:

1 X CBD Column Binding Buffer

Tris-HCl	20 mM Adjust pH to 8 using HCl	(2,42 g/L)
NaCl	500 mM	(29,22 g/L)
EDTA	1 mM	(0,292 g/L)
Tween-20	0.05% v/v	500 ul/ 1 L

Gently mix the beads solution and resuspend thoroughly.



Chitin Magnetic Beads - 5 ml

by New England Biolabs

Catalog #: E8036S

## 2 Aliquot 50 ul of beads solution into a 2 ml microfuge tube

## 3 Apply a magnet for 30 seconds to pull the beads to one side. Remove the excess liquid

03s

- 4 Wash the beads 3 times with 500 ul of 1x CBD Column Binding Buffer. Apply magnet and remove the liquid.
- 5 Load 200-500 ul of cell supernatant with the beads.  
Note: when making dilutions of the supernatant, dilute in CBD Column Binding Buffer.
- 6 Mix thoroughly and incubate for 1 hour at 4°C at constant agitation. 1h
- 7 Apply magnet for 30 seconds and remove supernatant with a pipette. Save the supernatant. 03s
- 8 Wash beads 3 times with 500 ul of 1x CBD Column Binding Buffer. Apply magnet for 30 seconds and pipette off liquid. Save the wash fractions. 03s
- 9 To remove proteins from the beads, boil beads for 5 minutes in 50 ul SDS-Page sample buffer at 95 °C 5m
- 10 Run fractions on a protein gel



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited