



Mar 28,  
2019

Working

## ROS staining for Arabidopsis (green fluorescent stain)

Magdalena Julkowska<sup>1</sup>

<sup>1</sup>King Abdullah University of Science and Technology

[dx.doi.org/10.17504/protocols.io.zksf4we](https://doi.org/10.17504/protocols.io.zksf4we)



Magdalena Julkowska  
King Abdullah University of Science and Technology



### ABSTRACT

This protocol is adapted from two publications describing fluorescent staining for ROS - [Shin et al., 2005](#) and [Orman-Ligeza et al. 2016](#) (citing Shin)

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS

NAME	CATALOG #	VENDOR
Murashige & Skoog medium including B5 vitamins	M0231	Duchefa Biochemie
Costar® 6 Well Clear TC-Treated Multiple Well Plates, Individually Wrapped, Sterile	3516	Westnet
CM-H2DCFDA	C6827	Thermo Fisher Scientific

### BEFORE STARTING

prepare sterile MS medium:

- 1/2 MS
- 0.5% sucrose
- 0.1% MES buffer
- pH 5.8 with KOH

Autoclave

prepare CM-H2DCFDA according to the supplier instructions

- 1 Sterilize the seedlings, germinate and transfer to for e.g. salt stress medium as described in the "Quantification of salt-induced changes in Root System Architecture in Arabidopsis" protocol - [dx.doi.org/10.17504/protocols.io.zkqf4vw](https://doi.org/10.17504/protocols.io.zkqf4vw)
- 2 Transfer the seedlings from agar plates into nutrient solution that is the same as agar plate composition (but without Dashin agar). You can use 24/12/6 well-plates - this will depend on your final seedlings size
- 3 After transferring the seedlings to liquid medium add DCF-DA so that the final concentration of the DCF-DA in the incubation medium is 50 uM
- 4 Wrap the well plate in aluminium foil (DCF-DA is light sensitive) and let the samples incubate for 30 min with gentle shaking.
- 5 Wash the seedlings with medium without DCF-DA

- 6 Image the stain using fluorescence microscopy:
- 460–500 nm bandpass excitation
  - 510– 560 nm bandpass emission



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited