

# Vacuole staining in *Phaeodactylum tricornutum* using RatioWorks™ PDMPO Version 2

Jernej Turnsek

## Abstract

Simple and rapid vacuole staining protocol for *Phaeodactylum tricornutum*. Note that RatioWorks™ PDMPO and LysoSensor™ Yellow/Blue DND-160 are the same thing and can thus be used interchangeably. RatioWorks™ PDMPO was used in this protocol.

**Citation:** Jernej Turnsek Vacuole staining in *Phaeodactylum tricornutum* using RatioWorks™ PDMPO. **protocols.io**  
dx.doi.org/10.17504/protocols.io.nbxdapn

**Published:** 18 Feb 2018

## Guidelines

## Materials

LysoSensor™ Yellow/Blue DND-160 L7545 by [Thermo Fisher Scientific](#)

RatioWorks™ PDMPO [21204](#) by [AAT Bioquest](#)

## Protocol

Prepare 5.46 mM RatioWorks™ PDMPO stock in 100% DMSO.

### Step 1.

Dilute 10x with 100% DMSO.

### Step 2.

Transfer 500 µL *Phaeodactylum tricornutum* culture to 1.5 mL epptube.

### Step 3.

Add 1 µL 10x RatioWorks™ PDMPO dilution to cells for ~1 µM final concentration.

### Step 4.

Incubate 10 min at room temperature.

**Step 5.**

Visualize with 405 nm laser line and emission window set to 525-575 nm.

**Step 6.**

Note: Leica TCS SP5 confocal microscope was used for images below.

📈 EXPECTED RESULTS

