

# QIAgen DNeasy Blood & Tissue kit, cultured cells

#### **Harold Bien**

### **Abstract**

QIAgen's DNeasy Blood & Tissue kit (cat#69504 or 69506) for cultured cells

Citation: Harold Bien QIAgen DNeasy Blood & Tissue kit, cultured cells. protocols.io

dx.doi.org/10.17504/protocols.io.dsz6f5

Published: 10 Sep 2015

### **Before start**

## Preparation of buffer AW1 and AW2

Add ethanol to AW1 and AW2 as directed on bottle. Stable for one year.

If buffer AL has precipitates, warm to 56C until precipitates dissolved

Preheat thermoixer to 56C

## **Equipment required:**

Pipettes and pipette tips Vortexer Microcentrifuge tubes, 1.5-1.7mL Thermomixer for heating at 56C

### **Materials**

✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users QIAgen DNeasy Blood and Tissue Kit, 50 rxn 69504 by Qiagen Ethanol, pure 4455 by Omnipure Filter Company Buffer AL, Lysis buffer 19076 by Qiagen Buffer AW1, Wash buffer (1), concentrate 19081 by Qiagen Buffer AW2, Wash buffer (2), concentrate 19072 by Qiagen Buffer AE, Elution buffer, 240mL 19077 by Qiagen

### **Protocol**

### Step 1.

Centrifuge a maximum of 5 million cells for 5 minutes at 300g

© DURATION

00:05:00

### Step 2.

Resuspend pellet in 200µL of PBS

**■** AMOUNT

200 μl Additional info:



✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users

### Step 3.

Add 20µL proteinase K.

**■** AMOUNT

20 μl Additional info:



Proteinase K, 2mL 19131 by Qiagen

### Step 4.

Add 200µL Buffer AL and mix thoroughly by vortexing

**■** AMOUNT

200 µl Additional info:



Buffer AL, Lysis buffer 19076 by Qiagen

### Step 5.

Add 200µL pure ethanol and mix by vortexing

**■** AMOUNT

200 μl Additional info:



Ethanol, pure 4455 by Omnipure Filter Company

#### Step 6.

Pipet mixture into DNeasy Mini spin column placed on a 2mL collection tube (supplied with kit)

### Step 7.

Centrifuge column at ≥6,000g for 1 minute

© DURATION 00:01:00

# Step 8.

Discard collection tube and flow through

#### Step 9.

Place spin column in new 2mL collection tube

### Step 10.

Add 500µL buffer AW1

AMOUNT

500 µl Additional info:



Buffer AW1, Wash buffer (1), concentrate 19081 by Qiagen

### **Step 11.**

Centrifuge for 1 minute at≥6,000g

**O DURATION** 

00:01:00

### **Step 12.**

Discard collection tube and flow through

# **Step 13.**

Add 500µL of Buffer AW2

**■** AMOUNT

500 μl Additional info:



Buffer AW2, Wash buffer (2), concentrate 19072 by Qiagen

### **Step 14.**

Centrifuge at 20,000g for 3 minutes

**O** DURATION

00:30:00

# Step 15.

Discard collection tube and flow through

### Step 16.

Transfer spin column to new (final and labeled) collection tube [optionally sterilized]

### Step 17.

Elute DNA by adding 200µL of buffer AE to the center of the spin column

**■** AMOUNT

200 µl Additional info:



**REAGENTS** 

Buffer AE, Elution buffer, 240mL 19077 by Qiagen

### **Step 18.**

Incubate for 1 minute at room temperature

© DURATION

00:01:00

### Step 19.

Centrifuge at ≥6,000g for 1 minute

**O** DURATION

00:01:00

### Step 20.

[Optiona]: Repeat step 19 for increased DNA yield