



Mar 25,  
2019

Working

## Polymerase Chain Reaction

Version 1

Lis Rocha<sup>1</sup>

<sup>1</sup>Universidade Federal de Viçosa

[dx.doi.org/10.17504/protocols.io.zhff33n](https://doi.org/10.17504/protocols.io.zhff33n)



Lis Rocha  
Universidade Federal de Viçosa



### ABSTRACT

This was the PCR protocol used for each set of primers in the study Comparative genomics of *Staphylococcus aureus* associated with subclinical and clinical bovine mastitis (Rocha et al., 2019)

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS TEXT

50 ng of total DNA, 1U of Taq DNA polymerase Cellco Biotec, 0.2  $\mu$ M of each primer, 0.2 mM deoxynucleotide triphosphate mixture, 1X reaction buffer containing 2.0 mM MgCl<sub>2</sub>, extra 1.0 mM MgCl<sub>2</sub>, and Milli-Q water to increase the reaction volume to a final volume of 25  $\mu$ L.

The extra 1 mM MgCl<sub>2</sub> was excluded from the PCR reaction that contained the primers LipoP-F-CS/LipoP-R-C.

**Table 1 - Primer Sequences for primers used in this Protocol**

cl3309subF	TGTTGTAGGAGGAACAATCC
cl3309subR	TTCTAATGTCAGCAACATGC
cl3309cliF	GCTATTCCTAGATGCACT
cl3309cliR	TTTAAAGTATGACATGAATG
cl3316F	ACGCAAAACCCCTTACTAGT
cl3316R	GCAACAACCTAGTAGGAGTGA
LipoP-F-CS	GYTTTGCGAAAACGTTAGAYATGTA
LipoP-R-C	TGCCTTCATCATTATTGGACCAATC
LipoP-F-CS	GYTTTGCGAAAACGTTAGAYATGTA
LipoP-R-CS	GGTAAAYTCAATGTCTTATRTCC

#### primers cl3309sub F/R

- 1 Initial denaturation: 95.0 °C for 5 min;
- 2 35 cycles of denaturation at 95.0 °C for 45 s,
- 3 Annealing: 55 °C for 45s

4 Extension: 72 °C for 45 s

5 final extension at 72.0 °C for 10 min

primers cl3316F/R

6 initial denaturation: 95.0 °C for 5 min;

7 35 cycles of denaturation at 95.0 °C for 45 s,

8 Annealing: 55 °C for 45 s

9 Extension: 72 °C for 45 s

10 final extension at 72.0 °C for 10 min.

primers cl3700 - LipoP FCS/RC

11 initial denaturation: 95.0 °C for 5 min;

12 35 cycles of denaturation at 95.0 °C for 45 s,

13 Annealing: 54 °C for 45 s

14 Extension: 72 °C for 45 s

15 final extension at 72.0 °C for 10 min.

cl33009cli F/R

16 initial denaturation: 95.0 °C for 5 min;

17 35 cycles of denaturation at 95.0 °C for 45 s,

18 Annealing: 45 °C for 45 s

19 Extension: 72 °C for 30 s

20 final extension at 72.0 °C for 10 min.

primers cl3700 - LipoP FCS/RCS

21 initial denaturation: 95.0 °C for 5 min;

22 35 cycles of denaturation at 95.0 °C for 45 s,

23 Annealing: 50 °C for 45 s

24 Extension: 72 °C for 1min

25 final extension at 72.0 °C for 10 min.

Analyzing the amplified fragments

26 Analyze the amplicons by electrophoresis in 1X Tris-acetate-EDTA on a 1.0% agarose gel and visualize imagen under UV light after staining with 2 mg.ml-1 ethidium bromide.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited