

Sep 25, 2019

Field sampling of root-associated microbes for DNA/RNA extraction V.1

Roey Angel¹¹Soil and Water Research Infrastructure

In Development

dx.doi.org/10.17504/protocols.io.qprdv6



Roey RA. Angel

Soil and Water Research Infrastructure



ABSTRACT

This protocol describes a procedure for sampling plant roots in the field for future DNA and RNA extraction for microbiome analysis. The protocol is deliberately designed to be simple and requires no electronic equipment. Root samples are preserved in LifeGuard Soil Preservation Solution for protecting against nucleic acid degradation.

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Micro-spatula set	AT16.1	Carl Roth
LifeGuard Soil Preservation Solution	12868-100	Qiagen
Scissors	HCT7.1	Carl Roth
Technical-grade ethanol (70%)	T913.1	Carl Roth
Paper towels	Y03.1	Carl Roth
Microcentrifuge tubes 2 ml	CK06.1	Carl Roth
Garden trowel	View	Amazon
Disposable pasteur pipettes	EA61.1	Carl Roth
Tweezers set	PX40.1	Carl Roth
Cooling box	AA46.1	Carl Roth
Cooling packs	E447.1	Carl Roth

STEPS MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Garden trowel	View	Amazon
Micro-spatula set	AT16.1	Carl Roth
Scissors	HCT7.1	Carl Roth
Paper towels	Y03.1	Carl Roth
Scissors	HCT7.1	Carl Roth
Paper towels	Y03.1	Carl Roth
Tweezers set	PX40.1	Carl Roth
Microcentrifuge tubes 2 ml	CK06.1	Carl Roth
Micro-spatula set	AT16.1	Carl Roth
LifeGuard Soil Preservation Solution	12868-100	Qiagen

NAME ▾	CATALOG # ▾	VENDOR ▾
Disposable pasteur pipettes	EA61.1	Carl Roth
Cooling box	AA46.1	Carl Roth
Cooling packs	E447.1	Carl Roth

BEFORE STARTING

Clean spatulas using 70% ethanol

G
M
T

<input type="checkbox"/>	Detect language	English					
--------------------------	-----------------	---------	--	--	--	--	--

Text-to-speech function is limited to 200 characters

<input checked="" type="checkbox"/>		Options : History : Feedback : Donate	Close
-------------------------------------	--	--	-------

- 1 **Sample a triplicate of plant individuals spaced out a few metres apart from each other** Make sure you are reall sampling individual plants and not offshoots of the same plant


- 2 **Using a garden trowel, carefully dig out the plants while keeping the root system intact** (as much as possible, of course)

G
M
T

<input type="checkbox"/>	Detect language	English					
--------------------------	-----------------	---------	--	--	--	--	--

Text-to-speech function is limited to 200 characters

<input checked="" type="checkbox"/>		Options : History : Feedback : Donate	Close
-------------------------------------	--	--	-------



Garden trowel
by Amazon
[View](#)

- 3 **While holding the plant by the shoot, shake the root system hard enough so that all loose soil is removed from it.** Take care to damage the plant as little as possible You can use a spatula to remove large soil aggregates that are attached to the roots



Micro-spatula set
by Carl Roth
Catalog #: [AT16.1](#)

- 4 **From the remaining root system (plus soil particles plus attached to the roots), trim a 'representative sample' of roots using scissors or scalpel** It is usually best to trim the roots onto a piece of paper towel



Scissors

by Carl Roth

Catalog #: [HCT7.1](#)



Paper towels

by Carl Roth

Catalog #: [Y03.1](#)

- 5 **Cut the trimmed out roots a little so that they fit into a 2.0 ml tube**



Scissors

by Carl Roth

Catalog #: [HCT7.1](#)



Paper towels

by Carl Roth

Catalog #: [Y03.1](#)

- 6 **Place about 2-3 g of that cut out sample into a 2.0 ml tube**



Tweezers set

by Carl Roth

Catalog #: [PX40.1](#)



Microcentrifuge tubes 2 ml

by Carl Roth

Catalog #: [CK06.1](#)

2 g



The root tissue should make up at least half or more of the mass, while the remaining attached soil should make up the rest

7 **Press the sample a little into the bottom of the tube to decrease its volume**



Micro-spatula set

by [Carl Roth](#)

Catalog #: [AT16.1](#)



Make sure the roots do not take up more than 1/2-2/3 of the volume

It is of course possible to split each sample into several separate tubes, depending on the specific type of roots, and submerge each with LifeGuard solution

8 **Add as much LifeGuard solution so that the sample is submerged in about twice of its volume (about 1.0 – 1.5 ml)** Best is to use a disposable Pasteur pipette for dispensing the solution



LifeGuard Soil Preservation Solution

by [Qiagen](#)

Catalog #: [12868-100](#)



Disposable pasteur pipettes

by [Carl Roth](#)

Catalog #: [EA61.1](#)

1.5 ml

9 **Place the tubes in cooling (around 4 °C) and keep them cooled until you reach the lab. The solution will protect nucleic acids even at room temperature for several days, but cooling is preferred**



Cooling box

by [Carl Roth](#)

Catalog #: [AA46.1](#)



Cooling packs

by [Carl Roth](#)

Catalog #: [E447.1](#)

4 °C

10 **In the lab, store the samples in a freezer (-20 – -80 °C)**

-20 °C or -80 °C



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited