

Milliplex Cytokine/Chemokine 7-plex MAG

Ligia Pinto, Troy Kem

Abstract

Milliplex Cytokine/Chemokine 7-plex manufacturer's protocol

Citation: Ligia Pinto, Troy Kem Milliplex Cytokine/Chemokine 7-plex MAG. **protocols.io**

dx.doi.org/10.17504/protocols.io.hvgb63w

Published: 01 Aug 2017

Protocol

Step 1.

PREPARATION OF SAMPLES/REAGENTS FOR IMMUNOASSAY

Step 2.

Preparation of Serum/Plasma Thaw Time: Thaw the samples completely on ice, mix well by shaking on plate shaker for 1 min. at RT (20-25°C) and centrifuge (1,700 xg, 10 minutes, 4°C) prior to use in the assay to remove particulates.

Step 3.

Preparation of Antibody-Immobilized Beads Sonicate each individual antibody-bead vial for 30 seconds; vortex for 1 minute. Add 60 µL from each antibody-bead vial to the Mixing Bottle and bring final volume to 3.0 mL with Bead Diluent. Vortex the mixed beads well. Unused portion may be stored at 2-8°C for up to one month. (Note: Due to the composition of magnetic beads, you may notice a slight color in the bead solution. This does not affect the performance of the beads or the kit). Example 1: When using 7 cytokine antibody-immobilized beads, add 60 µL from each of the 7 bead sets to the Mixing Bottle. Then add 2.58 mL Bead Diluent.

Step 4.

Preparation of Quality Controls Reconstitution Time: Before use, reconstitute Quality Control 1 and Quality Control 2 with 250 µL deionized water. Invert the vial several times to mix and vortex. Allow the vial to sit for 5-10 minutes and then transfer the controls to appropriately labeled polypropylene microfuge tubes. Unused portion may be stored at £ -20°C for up to one month.

Step 5.

Preparation of Wash Buffer Bring the 10X Wash Buffer to room temperature and mix to bring all salts into solution. Dilute 60 mL of 10X Wash Buffer with 540 mL deionized water. Store unused portion at 2-8°C for up to one month.

Step 6.

Preparation of Serum Matrix Reconstitution Time: Add 1.0 mL deionized water to the bottle containing lyophilized Serum Matrix. Mix well. Allow at least 10 minutes for complete reconstitution. Leftover reconstituted Serum Matrix should be stored at £ -20°C for up to one month.

Step 7.

Preparation of Human Cytokine Standard Reconstitution Time: 1.) Prior to use, reconstitute the Human Cytokine Panel III standard with 250 µL deionized water to give STD7. Invert the vial several

times to mix. Vortex the vial for 10 seconds. Allow the vial to sit for 5-10 minutes and then transfer the standard to an appropriately labeled polypropylene microfuge tube. This standard will be termed STD7; the unused portion may be stored at -20°C for up to one month.

Step 8.

) Preparation of Working Standards. Label six polypropylene microfuge tubes STD6, STD5, STD4, STD3, STD2, and STD1. Add 150 μL of Assay Buffer to each of the six tubes. -Prepare serial dilutions by adding 50 μL of STD7 reconstituted standard to the STD6 tube, mix well and transfer 50 μL of the STD6 standard to the STD5 tube, mix well and transfer 50 μL of the STD5 standard to the STD4 tube, mix well and transfer 50 μL of the STD4 standard to STD3 tube, mix well and transfer 50 μL of the STD3 standard to the STD2 tube, mix well and transfer 50 μL of the STD2 standard to the STD1 tube and mix well. The 0 pg/mL standard (Background) will be Assay Buffer. Standard Volume of Deionized Water to Add (mL) Volume of Standard to Add Original (STD7) 2500 Standard Concentration (pg/mL) Volume of Assay Buffer to Add (mL) Volume of Standard to Add STD6 150 50 mL of STD7 STD5 150 50 mL of STD6 STD4 150 50 mL of STD5 STD3 150 50 mL of STD4 STD2 150 50 mL of STD3 STD1 150 50 mL of STD2 Reconstituted STD6 STD5 STD4 STD3 STD2 STD1 Standard (STD7) Standard Tube # CXCL11/I-TAC (pg/mL) CXCL6/GCP2, CCL19/MIP3, CCL20/MIP3 (pg/mL) XCL1/Lymphotoxin, IL-11 (pg/mL) CXCL9/MIG (pg/mL) M-CSF (pg/mL) IL-29 (pg/mL) 10.49

Step 9.

4

Step 10.

9

Step 11.

2

Step 12.

4

Step 13.

82

Step 14.

Step 15.

8

Step 16.

5

Step 17.

8

Step 18.

7

Step 19.

33

Step 20.

8

Step 21.

1

Step 22.

1

Step 23.

3390.6

Step 24.

34

Step 25.

3

Step 26.

3

Step 27.

5

Step 28.

31,562.53,12551256251,2503,1256,25012,50065002,5005,00012,50025,00050,00072,00010,00020,00050,000100,000200,000IMMUNOASSAY PROCEDURE□ Allow all reagents to warm to room temperature (20-25°C) before use in the assay.□ Run the standards, controls, and samples in duplicate.

Step 29.

Prewet the plate by pipetting 200 µL of Assay Buffer into each well of the MAG Plate. Seal and mix on a plate shaker for 10 minutes at room temperature (20-25°C).

Step 30.

Decant Assay Buffer and remove residual amount from all wells by inverting the plate and tapping it smartly onto absorbent towels several times.10) Add 25 µL of each Standard or Control into the appropriate wells. Add 25 µL Assay Buffer to the 0 pg/mL standard (Background).

Step 31.

Add 25 µL of Assay Buffer to the sample wells.

Step 32.

Add 25 µL of the Serum Matrix solution to the background, appropriate standards, and control wells.

Step 33.

Add 25 µL of Sample into the appropriate wells.

Step 34.

Vortex Mixing Bottle and add 25 µL of the mixed Beads to each well. (Note: During addition of Beads, shake bead bottle intermittently to avoid settling. Due to the composition of magnetic beads, you may notice a slight color in the bead solution. This does not affect the performance of the beads or the kit).

Step 35.

Seal the plate with a plate sealer, cover it with the lid. Wrap a rubber band around the plate, lid and, shaker platform and incubate with agitation on a plate shaker 2 hours at room temperature (20-25°C).

Step 36.

Gently remove fluid by aspiration. Do not invert plate.

Step 37.

Wash plate 2 times with 200 µL/well of Wash Buffer, removing Wash Buffer by aspiration between each wash. To avoid washing/aspiration related bead loss, allow approximately 60 seconds between dispensing of the Wash Buffer and subsequent aspiration.

Step 38.

Add 25 µL of Detection Antibodies into each well (Note: Allow the Detection Antibodies to warm to room temperature prior to addition.)

Step 39.

Seal, cover with lid, and incubate with agitation on a plate shaker for 1 hour at room temperature (20-25°C). DO NOT ASPIRATE AFTER INCUBATION.20) Add 25 µL Streptavidin-Phycoerythrin to each well containing the 25 µL of Detection Antibodies.

Step 40.

Seal, cover with lid and incubate with agitation on a plate shaker for 30 minutes at room temperature (20-25°C).

Step 41.

Gently remove all contents by aspiration. Do not invert plate.

Step 42.

Wash plate 2 times with 200 µL/well of Wash Buffer, removing Wash Buffer by aspiration between each wash. To avoid washing/aspiration related bead loss, allow approximately 60 seconds between dispensing of the Wash Buffer and subsequent aspiration.

Step 43.

Add 150 µL of Sheath Fluid to all wells. Resuspend the beads on a plate shaker for 5 minutes.

Step 44.

Run plate on Luminex 100™ IS.

Step 45.

Save and analyze the data using Bio-Plex Manager software.EQUIPMENT SETTINGSEvents: 50, per bead region Sample Size: 100 µLGate Settings 5000 to 25,000Time Out 60 secondsQUALITY CONTROLThe ranges for each analyte in Quality Control 1 and 2 are provided on the card insert or can be located at the MILLIPORE website www.millipore.com/techlibrary/index.do using the catalog number as the keyword.NOTES:Legend: A- Hemolyzed sample; B- Lipemic sample; C- Sample missing; D- Clogged Filter well; E- Low bead count; F- >30% Bead aggregation; G- Instrument repeatProcedure-Quick Reference

Milliplex Cytokine/Chemokine 7-plex manufacturer's protocol

Step 46.