



Apr 09, 2019

In devel.

RNA Isolation for Tissue using TRIzol

Version 2

Sze-Xian Lim¹¹Duke University[dx.doi.org/10.17504/protocols.io.zvwf67e](https://doi.org/10.17504/protocols.io.zvwf67e)

Sze-Xian Lim

- 1 Add **1 ml TRIzol** of TRIzol per **30 mg tissue** of tissue and homogenize using handheld homogenizer.
- 2 Incubate at **25 °C** for **00:05:00** to allow nucleoprotein complexes to dissociate.
- 3 Add **200 µl chloroform** of chloroform to each ml of TRIzol carefully and vortex to mix well.
- 4 Centrifuge at max speed for **00:15:00** at **4 °C**.
- 5 Carefully remove the top aqueous phase and transfer to a new Eppendorf tube. The interphase and bottom organic phase can be saved for DNA and protein respectively.
- 6 Add **500 µl 100% isopropanol** of 100% isopropanol to the aqueous phase, mix by inversion and incubate at **-20 °C** for a minimum of **02:00:00**.
- 7 Spin down at maximum speed for 30 mins to precipitate RNA.
- 8 Remove supernatant, and add **1 ml** of **[M]75 Volume Percent ethanol** ethanol to wash pellet.
- 9 Spin down at max speed for 15 minutes and remove supernatant.

10 Resuspend pellet in appropriate volume of nuclease free water.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited