

DCX Immunohistochemistry rotocol

Nara Gyzely de Moraes Magalhães , Cristovam Guerreiro Diniz and Cristovam Wanderley Picanço Diniz

Abstract

Protocol Immunohistochemistry free-floating sections with anti-doublecortin antibody for avian tissue.

Citation: Nara Gyzely de Moraes Magalhães , Cristovam Guerreiro Diniz and Cristovam Wanderley Picanço Diniz DCX Immunohistochemistry rotocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.h4ub8ww

Published: 19 May 2017

Protocol

Antigen Retrieval Method

Step 1.

1. To remove the excess paraformaldehyde fixative, wash with 0,1 M PBS the sections free-floating at agitation in room temperature for 3 min, 3 times.
2. Incubate with 12% Boric acid (pH =9,0, 70°C) in water bath for 1 hour. When the temperature reaches 50°C, using a brush put the sections at the recipient and wait the temperature reaches 70°C and keep at this temperature for 1 hour. Remove the recipients from water bath and wait until it is at room temperature.
3. Wash the sections with 0,1% PBS/T for 5 minutes at , 3 times.
4. Wash the sections with 0,1M PBS for 2 minutes (shaking), 3 times each.

Protein Blocking Step

Step 2.

1. Incubate with blocking buffer (Normal Horse Serum 10% in 0,3% PBST) for 12 hours at gentle agitation in refrigeration 4°C.

Primary Antibody

Step 3.

1. Remove the serum and incubate with Anti Doublecourtin antibody (Doublecortin C-18, sc-8066 Santa Cruz Biotechnology – 1:500), diluted in 0,3%PBS/T at gentle agitation, overnight at 4°C.
2. Wash the sections with 0,1M PBS/T 0,1% for 2 minutes (shaking), 3 times each.

Secondary Antibody

Step 4.

1. Incubate with secondary antibody (Biotinylated Horse Anti-Goat IgG Antibody, BA-9500, Vector

Laboratories) Diluted in 0,3% PBS/T at 1:400, during 1 hour at room temperature.

Blocking Step

Step 5.

1. Incubate with 0,3% hydrogen peroxide (diluted in 0,1 M PBS) during 15 minutes with light shaking.
2. Wash sections with PBS/T 0,1% during 2 minutes, 3 times (shaking).

ABC

Step 6.

1. Incubate in VECTASTAIN® ABC KIT solution (first 37,5 µl A + 37,5µl B with 1,88 ml 0,3%PBS/T for 30 minutes, after add 13,12 ml 0,3%PBS/T) during 1 hour at 4°C with light shaking.
2. Wash with 0,1% PBS /T for 5 minutes, 2 times (shaking).

DAB Visualization

Step 7.

1. GDN preparation
 - Firstly prepare the Solution A by mix 0,006g of Diaminobenzidine (DAB) with 5 ml of distilled water.
 - Secondly prepare the Solution B by mixing 0,250g of Nickel ammonium sulfate with 5 ml de Acetate Buffer pH 6.0
 - Thirdly mix Solution A and B adding ammonium chloride (0,004g) with 0.020g α-D-Glucose.

Leave the section in this mix during 5 minutes.
2. Incubate sections with solution GND and wait for 3 minutes, after add 0,007g of Glucose-oxidase for each 3ml of GND solution for revelation. Stop revelation when the goal contrast is achieved (use a low gain microscope).
3. Remove the GDN + Glucose oxidase and wash the sections using 0,1M PBS for 3 times (2 minutes each time) with light shaking.
4. Mount the sections in appropriate gelatinized microscope slides and dry at room temperature for 12 hours or more depending on the mounting medium of choice.
5. Dihydrate and add the cover slips.



REAGENTS

- Boric acid [View](#) by [P212121](#)
- Sodium acetate [View](#) by [P212121](#)
- Ammonium Chloride [View](#) by [P212121](#)
- ✓ Sodium Phosphate monobasic by Contributed by users
- Triton X-100 [T8787-50ML](#) by [Sigma Aldrich](#)
- Sodium phosphate dibasic [7558-79-4](#) by [Sigma Aldrich](#)
- Normal Horse Serum [S-2000](#) by [Vector Laboratories](#)

Doublecortin Antibody [sc-8066](#) by [Santa Cruz Biotechnology](#)

✓ Biotinylated Horse Anti-Goat IgG Antibody [BA-9500](#) by Contributed by users

VECTASTAIN Elite ABC HRP Kit (Peroxidase, Standard) [PK-6100](#) by [Vector Laboratories](#)

Ammonium nickel(II) sulfate hexahydrate [7785-20-8](#) by [Sigma Aldrich](#)

3,3'-Diaminobenzidine tetrahydrochloride [D5905](#) by [Sigma Aldrich](#)

α -D-Glucose [492-62-6](#) by [Sigma Aldrich](#)

⊕ NOTES

Nara Magalhaes 19 May 2017

After each step remove with a pipette the remaining solution from previous step, with careful not damage or lose sections.

Use only sterilized material to minimize risk of contamination between different antibodies or solutions.

At the end add chlorine to inactivate GND + DAB residual solution before disposal on appropriate container.