

# Ca<sup>2+</sup> imaging (salivary glands) from Yao CK et al. (2017)

Chi-Kuang Yao, Yu-Tzu Liu, I-Chi Lee, You-Tung Wang, Ping-Yen Wu

## Abstract

This protocol is from 'Flower Ca<sup>2+</sup> channel in CME and ADBE' of Yao CK et al.

Please see the manuscript for the full method details.

**Citation:** Chi-Kuang Yao, Yu-Tzu Liu, I-Chi Lee, You-Tung Wang, Ping-Yen Wu Ca<sup>2+</sup> imaging (salivary glands) from Yao CK et al. (2017). **protocols.io**

[dx.doi.org/10.17504/protocols.io.hjwb4pe](https://dx.doi.org/10.17504/protocols.io.hjwb4pe)

**Published:** 04 Apr 2017

## Before start

You'll need:

### 0 mM Ca<sup>2+</sup> hemolymph-like (HL)-3 solution:

- 70 mM NaCl
- 5 mM KCl
- 10 mM MgCl<sub>2</sub>
- 10 mM NaHCO<sub>3</sub>
- 5 mM trehalose
- 5 mM HEPES (pH 7.2)
- 115 mM sucrose

### 2 mM Ca<sup>2+</sup>/5 mM K<sup>+</sup>/7 mM glutamate solution:

- 70 mM NaCl
- 5 mM KCl
- 10 mM MgCl<sub>2</sub>
- 10 mM NaHCO<sub>3</sub>
- 5 mM trehalose
- 5 mM HEPES (pH 7.2)

- 115 mM sucrose
- 2 mM  $\text{CaCl}_2$
- 7 mM Monosodium glutamate

**4  $\mu\text{M}$  Fluo-4 AM (Invitrogen)/100 mM  $\text{Ca}^{2+}$ /HL-3 solution to be loaded with Fluo-4 AM dye**

## Protocol

### Ca<sup>2+</sup> imaging in salivary glands

#### Step 1.

Dissect the gland cells of the third instar larvae in 0 mM  $\text{Ca}^{2+}$  HL-3 solution and plate on poly-L-Lysine-coating coverslip.

### Ca<sup>2+</sup> imaging in salivary glands

#### Step 2.

Bath the cells in 4  $\mu\text{M}$  Fluo-4 AM (Invitrogen)/100 mM  $\text{Ca}^{2+}$ /HL-3 solution to be loaded with Fluo-4 AM dye.

### Ca<sup>2+</sup> imaging in salivary glands

#### Step 3.

Capture the images every 10 min during dye loading using MetaMorph software and ANDOR iXon 897 camera.

### La<sup>3+</sup> treatment

#### Step 4.

Incubate the cells in 4  $\mu\text{M}$  Fluo-4 AM (Invitrogen)/100 mM  $\text{Ca}^{2+}$ /100 mM  $\text{LaCl}_3$ /HL-3 solution.

### La<sup>3+</sup> treatment

#### Step 5.

Capture the images every 10 min during dye loading using MetaMorph software and ANDOR iXon 897 camera.

### Counting

#### Step 6.

Count the fluorescence intensity in the salivary gland cells and surrounding cover slips.

### Step 7.

Calculate final fluorescence value by subtracting fluorescence intensity in the gland cells from the dye fluorescence intensity in the coverslips.

#### 📌 NOTES

**Chi-Kuang Yao** 04 Apr 2017

The fluorescence intensity of one salivary gland was used for each data value. Image processing was achieved using Image J and LSM Zen.