

# Colony PCR on *V. natriegens*

Daniel Bauersachs

## Abstract

**Citation:** Daniel Bauersachs Colony PCR on *V. natriegens*. **protocols.io**

[dx.doi.org/10.17504/protocols.io.pukdnuw](https://doi.org/10.17504/protocols.io.pukdnuw)

**Published:** 30 Apr 2018

## Protocol

### Step 1.

|         |                    |
|---------|--------------------|
| 5 µl    | reaction buffer    |
| 0,5 µl  | dNTPs              |
| 0,5 µl  | primer forward     |
| 0,5 µl  | primer reverse     |
| 0,1 µl  | one-taq polymerase |
| 18,4 µl | A.dest.            |
| = 25 µl |                    |
|         |                    |

Now you have to pick one single colony with a pipette tip and put it into your reaction mix. After 15 seconds you can take the tip out and start the PCR.