



Jun 05,
2019

Thawing and Seeding Frozen Cells

Version 2 ▼

Kenneth Schackart¹

¹University of Arizona

Working

[dx.doi.org/10.17504/protocols.io.3rcgm2w](https://doi.org/10.17504/protocols.io.3rcgm2w)

 Kenneth Schackart 

ABSTRACT

How to thaw cells from the liquid nitrogen storage and seed into a tissue culture flask

GUIDELINES

- Gloves must be worn at all times.
- Perform all tasks within biosafety cabinet.
- Anything entering biosafety cabinet must be generously sprayed with 70% ethanol (even you).
- When finished, wipe biosafety cabinet with 70% ethanol, and UV for at least 15 minutes.

MATERIALS TEXT

(1) T-75 or T-25 flask per frozen cell vial (or more if plating at a lower density)


(1) 15 mL centrifuge tube per frozen cell vial

(1) 10 mL serological pipet tip per T-75 flask or (1) 5 mL serological pipet tip per T-25 flask

Warmed cell culture media

1000 µL filter pipette tips


Thaw Cells

- 1 Thaw cells by suspending cryotube in  **37 °C** water bath until completely thawed, but no longer than necessary

Transfer cell suspension


- 2 Within biosafety cabinet, transfer cell suspension to 15 mL centrifuge tube using 1000 µL pipette.

Dilute freezing medium


- 3 Add  **1 ml** warmed cell culture medium to cell suspension *dropwise*.



Adding the initial cell culture medium slowly helps prevent cell death caused by a rapid change in osmotic pressure.

- 4 Add an additional  3 ml warmed cell culture medium to cell suspension slowly.

Centrifuge cell suspension


- 5 Centrifuge the cell suspension at  1500 rpm for  00:03:00 .

Resuspend Cells

- 6 Remove bulk of supernatant with serological pipet, then remove remainder with 1000 µL pipette.



For small cell pellets, you are better off leaving a small amount of media than disturbing the cell pellet.

- 7 Add  1 ml warmed cell culture media to cell pellet.



Allowing the cell pellet rest in media for about 2 minutes will help with resuspension.

- 8 Gently pipette mix the cell pellet into the solution.

- 9 Add an additional  7 ml warmed cell culture media [ 3 ml for T-25].

Seed Cells

- 10 Using a serological pipet, transfer the cell suspension to the tissue culture flask.

Label Flask

- 11 Label the flask with:
- Cell line
 - Passage number
 - Date
 - Your initials

Incubate

- 12 Transfer flask to CO₂ incubator.

13 Don't forget to remove the vial you used from the frozen storage inventory.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited