

# Culture of unicellular protists in marine broth

Maria Rubio-Brotons

## Abstract

**Citation:** Maria Rubio-Brotons Culture of unicellular protists in marine broth. **protocols.io**

dx.doi.org/10.17504/protocols.io.ggxbtxn

**Published:** 22 Nov 2016

## Protocol

Suspend 37.4 g of Marine Broth powder in 1 L of distilled water

### Step 1.

Marine Broth (Ref. Marine Broth 2216, BD-Difco)

Heat with frequent agitation and boil for 1 minute

### Step 2.

Autoclave for 15 minutes at 121°C

### Step 3.

Allow to cool and filter through 0.22µm membrane filter

### Step 4.

To perform inside the hood

Store at 4°C

### Step 5.