

# Low Gel Temperature Agarose (LGTA) Media Version 2

G Jason Smith, April Woods

## Abstract

Use for solid phase culture and chemical selection of *Pseudo-nitzschia spp* and other recalcitrant diatoms

**Citation:** G Jason Smith, April Woods Low Gel Temperature Agarose (LGTA) Media. **protocols.io**

[dx.doi.org/10.17504/protocols.io.gmnbu5e](https://doi.org/10.17504/protocols.io.gmnbu5e)

**Published:** 01 Dec 2016

## Guidelines

Reagents:

- 0.2 µm filtered seawater (FSW)
- L1 Medium Kit ([ncma.bigelow.org](http://ncma.bigelow.org))
- SeaPrep™ Agarose (Lonza Cat#: 50302, 25 gm)

## Protocol

PREPARATION: 1% (w/v) LGTA

### Step 1.

Add 100 mL FSW to autoclave compatible screw top bottle with stir bar.

PREPARATION: 1% (w/v) LGTA

### Step 2.

To stirring FSW add:

- 1.00 gm SeaPrep agarose to stirring FSW, *avoid clumping*
- 0.10 mL L1 - NaH<sub>2</sub>PO<sub>4</sub> stock (36.2 µM final, 1X L1 conc)
- 0.05 mL L1- NaNO<sub>3</sub> stock (441 µM final, 0.5X L1 conc)
- 0.20 mL L1 - Na<sub>2</sub>SiO<sub>3</sub> stock (212 µM final, 2X L1 conc)
- 0.10 mL L1 - Trace Element stock (1X L1 conc)

## NOTES

G Jason Smith 01 Dec 2016

Reagent sources listed under guidelines

PREPARATION: 1% (w/v) LGTA

**Step 3.**

Autoclave mixture to fully melt agarose.

PREPARATION: 1% (w/v) LGTA

**Step 4.**

Cool in water bath to 20° to 25°C with occasional stirring

PREPARATION: 1% (w/v) LGTA

**Step 5.**

Add 0.05 mL L1 – Vitamin stock (standard conc.), gently mix

🔗 **NOTES**

**G Jason Smith** 01 Dec 2016

**NOTE:** LGTA stocks can be held in molten form at  $\geq 20^{\circ}\text{C}$  or allowed to solidify at  $\leq 15^{\circ}\text{C}$

and remelted in microwave prior to cooling for use. *KEEP STERILE*

PLATING CULTURES:

**Step 6.**

Aliquot cooled, molten LGTA into sterile culture tubes (e.g. Falcon 2059).

- Volume should be sufficient to form thin layer (<5mm) in target plates or wells to enable in-gel imaging.
- *Recommend 0.2 mL per 24 well plate well, 2.5 mL for 30mm diameter petri plates, and 5 mL for 60 mm diameter plates.*

PLATING CULTURES:

**Step 7.**

Add selection reagents as needed.

PLATING CULTURES:

**Step 8.**

Add cell suspension **into** LGTA.

PLATING CULTURES:

**Step 9.**

Gently vortex to disperse cells.

#### PLATING CULTURES:

##### Step 10.

Transfer to culture plate avoiding introducing bubbles.

#### PLATING CULTURES:

##### Step 11.

Cool and incubate at 15°C to set agarose.

#### PLATING CULTURES:

##### Step 12.

Transfer to standard growth conditions as needed

- For species with high growth temperatures solidified LGTA can be transferred to appropriate incubation temperature after solidifying.
- For long term solid phase culture, the LGTA gel can be overlaid with L1 media.

#### PLATING CULTURES:

##### Step 13.

**NOTE:** Good gel strength is obtained at  $\geq 0.75\%$  (w/v) LGTA and optical clarity obtained at  $\leq 2\%$  LGTA).

We have had good luck with *Pseudo-nitzschia* growing in  $\geq 0.75\%$  to  $\leq 1\%$  LGTA.

Concentration of stock LGTA can be adjusted to accommodate larger culture volumes but a

1:1 (v/v) mixture of culture to LGTA should be considered the maximum mixing ratio for good gel formation.

