

# PCR of mouse PCSK9 from cDNA (cleavage template)

Noreen Wauford

## Abstract

**Citation:** Noreen Wauford PCR of mouse PCSK9 from cDNA (cleavage template). **protocols.io**

dx.doi.org/10.17504/protocols.io.e8xbhxn

**Published:** 28 Jun 2016

## Protocol

### Step 1.

## PCR mix

50 uL H<sub>2</sub>O

16 uL GC buffer

1.6 uL DNTPs

4 uL 10 uM M13 fwd (-41) primer

4 uL 10 uM M13 rev (-48) primer

1 uL cDNA

.8 uL Phusion

Divide into PCR tubes with 20 or 40 uL each

### Step 2.

## Thermocycler protocol

98\* 30 s

98\* 10 s

61\* 10 s

72\* 10 s

repeat steps 2-4 23 times

72\* 2 min

4\* hold

### **Step 3.**

## **Agarose gel**

Add 10 uL 6x loading buffer per 40 uL PCR reaction and vortex

Run on agarose gel with largest comb, loading 25 uL per well

Run @ 110 V for 30 min

### **Step 4.**

## **Cut gel**

Use a razor to carefully excise the bands (2 kb) under UV light.

### **Step 5.**

## **Extract from gel**

Load gel slices onto the Ultrafree DA columns and spin at 5,000 g for 10 min

### **Step 6.**

## **Lyophilize overnight**

### **Step 7.**

## **Add water**

Start with 20 uL, dilute (1.5 uL : 8.5 uL H<sub>2</sub>O), and run on PAGE to estimate concentration