

Extraction RNA of cells by Trizol

Maria Cármén Sales, Maryana Branquinho, Maysa Silva

Abstract



Citation: Maria Cármén Sales, Maryana Branquinho, Maysa Silva Extraction RNA of cells by Trizol. **protocols.io**
dx.doi.org/10.17504/protocols.io.ssdeea6

Published: 21 Aug 2018

Before start

Clean the work area with 70% alcohol. Use all filters and autoclaved tips. Refrigerate your centrifuge to 4C.

Materials

- ✓ RNase-free Water by Contributed by users
-  TRizol Reagent [15596026](#) by [Thermo Fisher Scientific](#)
- ✓ Isopropanol by Contributed by users
-  Chloroform 319988 by [Sigma](#)
- ✓ 0.5% Sodium dodecyl sulfate solution by Contributed by users
- ✓ Ethanol 75% by Contributed by users
- ✓ 0.1 mM EDTA - Ethylenediaminetetraacetic acid by Contributed by users

Protocol

Step 1.

Remove the culture medium from the plate and add 750µl of Trizol to 1×10^5 - 10^7 cells.

Homogenize with the pipette and transfer to a tube.

In this step, you can freeze this sample for up to 6 months in the -80C freezer or continue the protocol.

Step 2.

Incubate for 5 minutes to permit complete dissociation of the nucleoproteins complex.

Step 3.

Add 200 µl of chloroform per 1 mL of Trizol and incubate for 2 to 3 minutes.

Step 4.

Centrifuge the sample for 15 minutes at 12,000 xg at 4°C.

The mixture separates into a lower red phenol-chloroform, and interphase, and a colorless upper aqueous phase.

Step 5.

Transfer the aqueous (incolor) phase containing the RNA to a new tube by angling the tube at 45° and pipetting the solution out.

Step 6.

Add 500 µl of isopropanol per 1 mL of Trizol. Incubate for 10 minutes.

Step 7.

Centrifuge for 10 minutes at 12,000 × g at 4°C.

Step 8.

Discard the supernatant with a micropipetto and resuspend the pellet in 1 mL of 75% ethanol per 1 mL of Trizol.

Step 9.

Vortex the sample briefly and centrifuge for 5 minutes at 7500 × g at 4°C.

Step 10.

Discard the supernatant with a micropipettor. Vacuum or air dry the RNA pellet for 5–10 minutes.

Step 11.

Resuspend the pellet in 20–50 µL of RNase-free water, 0.1 mM EDTA, or 0.5% SDS solution by pipetting up and down.

Step 12.

Incubate in a water bath or heat block set at 55–60°C for 10–15 minutes.

Step 13.

Store in freezer -80°C until use or make the cDNA reaction then.

Warnings

The triazole reagent is toxic, gloves, lab coat, mask should be used. Manipulate the trizol in the exhaust hood.