

Phosphate-buffered Saline (PBS)

katy monteith

Abstract

This protocol describes how to make 1 L of PBS to a final concentration of 1x or 10x solution. The solution can be made RNase free with the addition of DEPC (diethyl pyrocarbonate) as in step 5.

Citation: katy monteith Phosphate-buffered Saline (PBS). **protocols.io**

dx.doi.org/10.17504/protocols.io.p4rdqv6

Published: 14 May 2018

Materials

Hydrochloric acid [H1758](#) by [Sigma Aldrich](#)

Sodium chloride [S7653](#) by [Sigma Aldrich](#)

Potassium chloride [P9333](#) by [Sigma Aldrich](#)

Disodium phosphate [S7907](#) by [Sigma Aldrich](#)

Monopotassium phosphate [P9791](#) by [Sigma Aldrich](#)

Diethyl pyrocarbonate [D5758](#) by [Sigma Aldrich](#)

Protocol

Step 1.

For a final concentration of 1x PBS weigh out the following amount of each reagent: NaCl- 8g, KCl- 0.2 g, Na₂HPO₄- 1.44 g, KH₂PO₄- 0.24 g. For a final concentration of 10x PBS weigh out the following amount of each reagent: NaCl- 80g, KCl- 2 g, Na₂HPO₄- 14.4 g, KH₂PO₄- 2.4 g.

Step 2.

Dissolve reagents in 800 mL of triple distilled H₂O using a magnetic stirrer and flea or by manual shaking.

Step 3.

Using a pH meter, adjust pH of solution to 7.4 (or 7.2 if required) by adding concentrated hydrochloric acid (HCl).

Step 4.

Add triple distilled H₂O to bring final volume of solution to 1 L

Step 5.

If require RNase free PBS, add 1 mL of DEPC (diethyl pyrocarbonate) per 1 L of PBS.

Step 6.

Autoclave.

Warnings

Read appropriate SDS for sodium chloride, potassium chloride, disodium phosphate, monopotassium phosphate, hydrochloric acid and diethyl pyrocarbonate before commencing work.

Sort By Relevance ▼