



Apr 23, 2019

In devel.

## Transfection of *Cryptocodinium cohnii* using labelled DNA

José A. Fernández Robledo<sup>1</sup><sup>1</sup>Bigelow Laboratory for Ocean Sciences

dx.doi.org/10.17504/protocols.io.z26f8he

Protist Research to Optimize Tools in Genetics (PROT-G)

José A. Fernández Robledo ⚡ 🌈 🦋

### ABSTRACT

José A. Fernández Robledo

Nastasia J. Freyria

Bigelow Laboratory for Ocean Sciences

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

In progress, unpublished

### Reagents

#### 1 Liquid medium

- L1 + 10 gr/L glucose NPN (<https://ncma.bigelow.org/ccmp316#.XLOpDZNKi2w>)
- L1 medium with nitrogen limited (without yeast extract and glucose) to obtain cyst from *C. cohnii* swimming cells

#### Plates

- Prepare Petri dishes containing L1 medium bacteriological agar (1.5%)
- Prepare PEG 8000 solution (0.4%) with L1 medium, sterilized by using a filter (0.22µm)
- Pour 1 ml of L1-PEG8000 solution onto dry agar plates; alternatively, pour 12.5 ml of L1-PEG8000 solution onto dry agar plates
- Let the PEG800 solution infuse for 24 h
- Store at 4°C until use

#### Plasmid vectors

- Propagate *pPmMOE[MOE]:GFP-11* (Cold et al., 2016; Fernández Robledo et al., 2008) in JM109 for miniprep



*pPmMOE[MOE]:GFP-11*

- Linearize 20 µg with *NotI* (e.g., New England Biolabs® inc., Ipswich, MA, USA), clean the restriction digestion (e.g., GenElute® PCR Clean-up kit)
- Store at -20°C until use

#### DNA probe construction and DNA labeling

- \*Use forward (5' -CCGCACATGTATGGTGAGCAAGGGCGAGGAGC- 3') the reverse (5' -CGTAGGACATGTCTGTACAGCTCGTCCATGCCG- 3') primers targeting *pPmMOE[MOE]:GFP-11* (Fernández Robledo et al., 2008) to amplify a 739 bp DNA
- Clean the amplicon as above and label 1 µg using Ulysis™ Alexa Fluor® 488 Nucleic Acid Labeling Kit (ThermoFisher Scientific) following the manufacturer's instructions
- Store labelled DNA were stored at -20°C until use

\*Alternatively use any other set of primers/target DNA to generate the amplicon for labeling

#### References

Cold, E.R., Freyria, N.J., Martínez Martínez, J., Fernández Robledo, J.A., 2016. An agar-based method for plating marine protozoan parasites of the genus *Perkinsus*. PLoS One 11, e0155015.

Fernández Robledo, J.A., Lin, Z., Vasta, G.R., 2008. Transfection of the protozoan parasite *Perkinsus marinus*. Mol Biochem Parasitol 157, 44-53.

#### Cell Preparation

### 2 Liquid culture

- Incubate at 24°C in the dark with 125 rpm in a rotatory shaking incubator; the doubling time of *C. cohnii* cells under replete nitrogen

conditions is 24-48 hours

- Harvest cells at log phase
- Centrifuge cells at 2,500 g for 10 min at room temperature
- Discard supernatant and maintain cells on ice until use

### Spheroplasts induction

*Cryptocodinium cohnii* spheroplasts are derived from both swimming cells and cyst cells (Kwok et al., 2007; Pozdnyakov et al., 2014).

- Prepare PEG8000 (20%) (wt/vol) in L1 medium
- Resuspend *C. cohnii* swimming cells with 1 ml of 20% L1-PEG8000 solution and vortex for 10 min
- Spread 1 ml of *C. cohnii* swimming cells or cysts on Petri dishes containing L1 medium bacteriological agar and L1-PEG8000 solution (0.4%)
- Incubate for 2 days at 28°C
- Cover the colonies on the plate with L1 medium, pour off the medium
- Pour new L1 medium and retain the second elutant containing the spheroplasts
- Count spheroplasts using a haemocytometer (x3)

### References

Kwok, A.C., Mak, C.C., Wong, F.T., Wong, J.T., 2007. Novel method for preparing spheroplasts from cells with an internal cellulosic cell wall. *Eukaryot Cell* 6, 563-567.

Pozdnyakov, I., Matantseva, O., Negulyaev, Y., Skarlato, S., 2014. Obtaining spheroplasts of armored dinoflagellates and first single-channel recordings of their ion channels using patch-clamping. *Marine drugs* 12, 4743-4755.

### Electroporation

- 3
- Collect *Cryptocodinium cohnii* cells in log phase or spheroplasts ( $1.0 \times 10^6$  cells/mL)
  - Electroporated with 5  $\mu$ g (2.5  $\mu$ g supercoiled, 2.5  $\mu$ g *NotI* linearized) of pPmMOE[MOE]GFP-11 using Amaxa®Cell Line Optimization Nucleofector™ solution V and program X-001 in a Nucleofector™ (Lonza)
  - Recover cells from the cuvette and incubate in L1 medium
  - As control for plasmid and electroporation, electroporate *P. marinus* PRA240 with pPmMOE[MOE]GFP-11 as reported elsewhere (Fernández Robledo et al., 2008)
  - Monitor cells for green fluorescence using standard FITC excitation/emission filters (488/507nm) under a transmitted-light fluorescence and/or confocal microscope


### Lipofection

- 4
- Lipofectamine® 3000 Reagent kit (ThermoFisher Scientific) in tissue culture plate 24-wells and 6-wells plates containing respectively 2-3 mL of L1 medium
  - Use 15  $\mu$ l of Lipofectamine and 1  $\mu$ g of labeled DNA in  $1.0 \times 10^6$  cells
  - Monitor cells for green fluorescence using standard FITC excitation/emission filters (488/507nm) under a transmitted-light fluorescence and/or confocal microscope

5



*Cryptocodinium cohnii* viability after electroporation





This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited