

TE Buffer

Ms Alex Aitken

Abstract

TE is a commonly used buffer solution in molecular biology, especially in procedures involving DNA or RNA. "TE" is derived from its components: Tris, a common pH buffer, and EDTA, a molecule that chelates cations like Mg^{2+} . The purpose of TE buffer is to solubilize DNA or RNA, while protecting it from degradation: 10 mM Tris, bring to pH 8.0 with HCl, 1 mM EDTA

Citation: Ms Alex Aitken TE Buffer. [protocols.io](https://doi.org/10.17504/protocols.io.c5uy6v)

[dx.doi.org/10.17504/protocols.io.c5uy6v](https://doi.org/10.17504/protocols.io.c5uy6v)

Published: 25 Jan 2016

Guidelines

| | |
|----------------|-----------|
| Component | For 500ml |
| 1M Tris pH 8* | 5ml |
| 0.5M EDTA pH 8 | 1ml |
| dH2O | 496ml |

*the pH is usually adjusted to RNA 8.0 for DNA and 7.5 forThe respective DNA and RNA nucleases are supposed to be less active at these pH values, but pH 8.0 can safely be used for storage of both DNA and RNA. EDTA further inactivates nucleases, by binding to metal cations required by these enzymes. Genomic and plasmid DNA can be stored in TE Buffer at 4°C (39.2°F) for short-term use, or -20°C (-4°F) to -80°C (-112°F) for long-term storage.

Protocol

Step 1.

Prepare 5ml of 1M Tris

📄 [AMOUNT](#)

61 g Additional info:

📄 [PROTOCOL](#)

. [1M Tris](#)

CONTACT: [VERVE Team](#)

Step 1.1.

60.57 g of Tris (hydroxymethyl) aminomethane in 0.5L Milli-Qwater

Step 1.2.

Bring pH to 8.0 using HCl

Making 0.5M EDTA

Step 2.

Prepare 1ml 0.5M EDTA

 [AMOUNT](#)

19 g Additional info:

 [PROTOCOL](#)

. [0.5M EDTA](#)

CONTACT: [VERVE Team](#)

Step 2.1.

Mix 18.6 g EDTA in 100ml Milli-Q water

 [AMOUNT](#)

19 g Additional info:

Step 2.2.

Bring pH to 8.0 using NaOH

 [NOTES](#)

VERVE Team 22 Jun 2015

EDTA will not be soluble until pH reaches 8.0–this will take time (hours)

Making 0.5M EDTA

Step 3.

Combine 5ml of 1M Tris and 1ml of 0.5M EDTA in 496ml dH₂O

Making 0.5M EDTA

Step 4.

Autoclave to sterilize.

 [NOTES](#)

VERVE Team 22 Jun 2015

It is best to remove a working aliquot and not repeatedly access the stock

Warnings

Repeated freeze-thaw cycles should be avoided