



Oct 19, 2019

## Hybridization of DNA oligos

[Laura Sánchez<sup>1</sup>](#)<sup>1</sup>AEGIS - Madrid iGEM Team 2019**1** *Works for me* [dx.doi.org/10.17504/protocols.io.8guhtww](https://doi.org/10.17504/protocols.io.8guhtww)**AEGIS - Madrid iGEM 2019****Manuela de las Casas**

Universidad Complutense de Madrid, AEGIS - Madrid iGEM 2019



### ABSTRACT

This protocol is used to hibridate two complementary DNA chains.

### MATERIALS TEXT

#### Equipment

- Thermocycler
- Vortex
- Centrifuge

#### Reagents

- 1.0 M NaCl
- Sterile H<sub>2</sub>O

#### DNA

- Promoters 5'
- Promoters 3'

### SAFETY WARNINGS

Lab coat and gloves should be weared throughout the whole experiment. All working surfaces must be clean and all the reactives should be treated following manufacturer's instructions.

#### Annealing of oligonucleotides

- 1 Add the required volume of H<sub>2</sub>O to the lyophilized oligonucleotides to obtain a concentration of 100 µM.  
Vortex both tubes for 30 s and incubate them at RT for 5 min to dissolve them.
- 2 Prepare the annealing mix by adding into a PCR tube:
  - > 45.5 µl of H<sub>2</sub>O
  - > 2.5 µl of 1.0 M NaCl
  - > 1 µl of oligo 5' (100 µM)
  - > 1 µl of the oligo 3' (100 µM)
- 3 Place the PCR tube with the mix in a thermocycler with the following annealing programme:
  - 5 min at 95°C
  - 1 min at 95°C
  - ramp down 1°C per cycle for 72 cycles
  - end by keeping the temperature at 10°C
- 4 Take 10 µl of the annealed oligonucleotides and dilute it with 90 µl of H<sub>2</sub>O to obtain a 0.2 µM concentration.

5 The annealed oligonucleotides stocks can be stored at  $-20^{\circ}\text{C}$  for future use.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited