

Fe-EDTA Stock solution for ESAW Media for Marine Phytoplankton

Ashley Humphrey

Abstract

Fe-EDTA Stock Solution for ESAW Media for Marine Phytoplankton

Citation: Ashley Humphrey Fe-EDTA Stock solution for ESAW Media for Marine Phytoplankton. **protocols.io**
dx.doi.org/10.17504/protocols.io.getbten

Published: 16 Nov 2016

Protocol

Fe-EDTA Stock solution, amounts per one Liter (1000 mL)

Step 1.

Na₂EDTA

 **AMOUNT**

4 g Additional info:

Fe-EDTA Stock solution, amounts per one Liter (1000 mL)

Step 2.

FeCl₃

 **AMOUNT**

3 g Additional info:

Fe-EDTA Stock solution, amounts per one Liter (1000 mL)

Step 3.

Bring solution to 1000 mL with MilliQ H₂O and filter sterilize with a 0.2µm filter. Store in the dark at 4°C