

# Generation of DNA fragments by sonication

William H. Wilson and Declan Schroeder

## Abstract

**Citation:** William H. Wilson and Declan Schroeder Generation of DNA fragments by sonication. **protocols.io**  
dx.doi.org/10.17504/protocols.io.dy57y5

**Published:** 18 Jan 2016

## Protocol

### Step 1.

Generation of DNA fragments by sonication is performed by placing a microcentrifuge tube containing the buffered DNA sample into an ice-water bath in a cup-horn sonicator.

### Step 2.

Sonication is conducted for a varying number of 10 s bursts using maximum output and continuous power.

### Step 3.

Exact conditions for sonication should be empirically determined for a given DNA sample before a preparative sonication is performed.

### Step 4.

Typically, 100 µg DNA in TE buffer is split into 10 aliquots of 35 µL.

### Step 5.

5 aliquots are subjected to sonication for increasing numbers of 10 s bursts.

### Step 6.

Aliquots from each time point are run on an agarose gel to determine optimal-sized DNA fragments (1–4 kb).

### Step 7.

Once optimal sonication conditions are determined, the remaining 5 aliquots (approximately 8 µg) are sonicated according to those predetermined conditions.

## NOTES

**Declan Schroeder** 12 Oct 2015

DNA can be blunt-ended and size-selected prior to downstream cloning.