

# NEBNext® Ultra™ DNA Library Prep Protocol for Illumina® With Size Selection (E7370)

#### **Isabel Gautreau**

## **Abstract**

This protocol is a group of the following 5 methods from the NEB website:

- 1. NEBNext End Prep
- 2. Adaptor Ligation
- 3. Size Selection of Adaptor-ligated DNA
- 4. PCR Amplification
- 5. Cleanup of PCR Amplification

Citation: Isabel Gautreau NEBNext® Ultra™ DNA Library Prep Protocol for Illumina® With Size Selection (E7370).

protocols.io

dx.doi.org/10.17504/protocols.io.cg5ty5

Published: 10 Feb 2015

## **Guidelines**

The NEBNext Ultra DNA Library Prep Kit for Illumina contains enzymes and buffers that are ideal to convert a small amount of DNA input into indexed libraries for next-generation sequencing on the Illumina platform (Illumina, Inc).

## **Adaptor Ligation:**

If DNA input is < 100 ng, dilute the NEBNext Adaptor for Illumina\* 1:10 in sterile water use immediately to a final concentration of 1.5  $\mu$ M.

#### **Size Selection:**

This protocol is suitable for input amounts of 50ng or higher. For lower input amounts we recommend that the protocol without size selection be followed.


This size selection is for libraries with 200 bp inserts only and is based on a starting volume of 100  $\mu$ l. For libraries with different size fragment inserts, refer to Table 1 for the appropriate volume of beads to be added.


Table 1: Recommended conditions for bead based size selection.

LIBRARY	APPROXIMATE INSERT SIZE	150 bp	200 bp	250 bp	300-400 bp	400-500 bp	500-700 bp
PARAMETERS	Total Library Size (insert + adaptor)	270 bp	320 bp	400 bp	400-500 bp	500-600 bp	600-800 bp
	1st Bead Selection	65	55	45	40	35	30
BE ADDED (μl)	2nd Bead Selection	25	25	25	20	15	15

Note: Be sure not to transfer any beads. Trace amounts of bead carry over may affect the optimal performance of the polymerase used in the NEBNext High-Fidelity 2X PCR Master Mix in the subsequent PCR step.

Figure 1: Examples of libraries prepared with human IMR-90 gDNA.





## **Materials**

NEBNext Ultra DNA Library Prep Kit for Illumina - 24 rxns <u>E7370S</u> by New England Biolabs

## **Protocol**

## **NEBNext End Prep**

#### Step 1.

Mix the following components in a sterile nuclease-free tube (Total volume  $65 \mu l$ ):

**₽** PROTOCOL

. NEBNext End Prep Mixture

**CONTACT:** Isabel Gautreau

**P** NOTES

Isabel Gautreau 20 Oct 2014

Starting Material: 5 ng-1 µg fragmented DNA.

Step 1.1.

End Prep Enzyme Mix 3.0 µl

#### **ANNOTATIONS**

## Peter Hoyt 12 Oct 2015

Will be optimizing for long read instruments.

## Step 1.2.

End Repair Reaction Buffer (10X) 6.5 μl

## **ANNOTATIONS**

## Peter Hoyt 12 Oct 2015

Might be in the wrong protocol as the Ultrall End Prep says use 7ul of reaction buffer.

## Step 1.3.

Fragmented DNA 55.5 µl

## **ANNOTATIONS**

## Peter Hoyt 12 Oct 2015

Initial experiment shooting for fragments greater than 20kb. Not sure why this protocol doesn't show the recommended 500pg - 1ug of fragmented DNA, but we will need approximately 4ug of HMW (High Molecular Weight) DNA to make of the molar difference in molecule ends.

## **NEBNext End Prep**

#### Step 2.

Mix by pipetting

## **NEBNext End Prep**

#### Step 3.

Quick spin to collect all liquid from the sides of the tube.

## **NEBNext End Prep**

#### Step 4.

Place in a thermocycler, with the heated lid on, and run the following program:

Time	Temperature			
30 minutes	20°C			
30 minutes 65°C				
Hold	4°C			

## **Adaptor Ligation**

#### Step 5.

Add the following components directly to the End Prep reaction mixture and mix well (Total volume  $83.5 \mu l$ ):

## **№** PROTOCOL

## . Adaptor Ligation Mixture (E7370)

**CONTACT: New England Biolabs** 

## NOTES

#### Isabel Gautreau 26 Jan 2015

The NEBNext adaptor is provided in NEBNext Singleplex (<u>NEB #E7350</u>) or Multiplex (<u>NEB #E7335</u>, <u>#E7500</u>) Oligos for Illumina.

## Step 5.1.

Blunt/TA Ligase Master Mix 15 μl

#### Step 5.2.

NEBNext Adaptor for Illumina 2.5 μl

#### NOTES

## New England Biolabs 29 Sep 2014

The NEBNext adaptor is provided in NEBNext Singleplex (<u>NEB #E7350</u>) or Multiplex (<u>NEB #E7335</u>, <u>#E7500</u>) Oligos for Illumina.

## **New England Biolabs** 20 Oct 2014

If DNA input is < 100 ng, dilute the NEBNext Adaptor for Illumina 1:10 in sterile water. Use immediately to a final concentration of 1.5  $\mu$ M.

#### **ANNOTATIONS**

## Peter Hoyt 12 Oct 2015

Not sure we need this step. Need to get back to this step to see if it helps create capped (ligated) ends.

## Step 5.3.

Ligation Enhancer 1 μl

## **Adaptor Ligation**

## Step 6.

Mix by pipetting

## **Adaptor Ligation**

## Step 7.

Quick spin to collect all liquid from the sides of the tube.

## **Adaptor Ligation**

#### Step 8.

Incubate at 20°C for 15 minutes in a thermal cycler.

#### **O DURATION**

00:15:00

## **Adaptor Ligation**

#### Step 9.

Add 3 µl of USER™ enzyme to the ligation mixture from step 8.



3 μl Additional info:



USER Enzyme - 50 units <u>M5505S</u> by <u>New England Biolabs</u>

## **Adaptor Ligation**

## Step 10.

Mix well and incubate at 37°C for 15 minutes.

**O DURATION** 

00:15:00

#### Size Selection

## **Step 11.**

Vortex AMPure XP beads to resuspend.

#### NOTES

## Isabel Gautreau 20 Oct 2014

The size selection protocol that follows is based on a starting volume of 100 µl.

Isabel Gautreau 20 Oct 2014

The following size selection protocol is for libraries with 200 bp inserts only. For libraries with different size fragment inserts, refer to Table 1 in the guidelines for the appropriate volume of beads to be added.

#### Size Selection

## **Step 12.**

Add 13.5 µl dH2O to the ligation reaction for a 100 µl total volume.

#### Size Selection

## **Step 13.**

Add 55 µl of resuspended AMPure XP beads to the 100 µl ligation reaction.

#### NOTES

#### Isabel Gautreau 20 Oct 2014

 $55~\mu l$  is the appropriate amount for libraries with 200 bp inserts only. For libraries with different size fragment inserts, refer to Table 1 in the guidelines for the appropriate volume of beads to be added.

#### Size Selection

#### **Step 14.**

Mix well by pipetting up and down at least 10 times.

#### Size Selection

#### Step 15.

Incubate for 5 minutes at room temperature.

**O DURATION** 

00:05:00

#### Size Selection

## **Step 16.**

Quickly spin the tube.

## Size Selection

#### Step 17.

Place the tube on an appropriate magnetic stand to separate the beads from the supernatant.

**O DURATION** 

00:05:00

#### Size Selection

## Step 18.

After the solution is clear (about 5 minutes), carefully transfer the supernatant containing your DNA to a new tube (**Caution: do not discard the supernatant**).

## Size Selection

## Step 19.

Discard the beads that contain the unwanted large fragments.

## Size Selection

#### Step 20.

Add 25 µl resuspended AMPure XP beads to the supernatant.

#### O NOTES

#### Isabel Gautreau 20 Oct 2014

 $25 \,\mu$ l is the appropriate amount for libraries with 200 bp inserts only. For libraries with different size fragment inserts, refer to Table 1 in the guidelines for the appropriate volume of beads to be added.

## Size Selection

#### Step 21.

Mix well and incubate for 5 minutes at room temperature.

© DURATION

00:05:00

#### Size Selection

## Step 22.

Quickly spin the tube.

#### Size Selection

#### Step 23.

Place it on an appropriate magnetic stand to separate the beads from the supernatant.

© DURATION

00:05:00

## Size Selection

#### Step 24.

After the solution is clear (about 5 minutes), carefully remove and discard the supernatant that contains unwanted DNA. Be careful not to disturb the beads that contain the desired DNA targets (**Caution: do not discard beads**).

#### Size Selection

#### **Step 25.**

Wash #1: Add 200 µl of 80% freshly prepared ethanol to the tube while in the magnetic stand.

## Size Selection

## Step 26.

Wash #1: Incubate at room temperature for 30 seconds

**O DURATION** 

00:00:30

## Size Selection

#### Step 27.

Wash #1: Carefully remove and discard the supernatant.

#### Size Selection

## **Step 28.**

Wash #2: Add 200 µl of 80% freshly prepared ethanol to the tube while in the magnetic stand.

## Size Selection

#### Step 29.

Wash #2: Incubate at room temperature for 30 seconds

## Size Selection

## Step 30.

Wash #2: Carefully remove and discard the supernatant.

#### Size Selection

#### **Step 31.**

Wash #3: Add 200 µl of 80% freshly prepared ethanol to the tube while in the magnetic stand.

## Size Selection

## Step 32.

Wash #3: Incubate at room temperature for 30 seconds

#### Size Selection

#### Step 33.

Wash #3: Carefully remove and discard the supernatant.

#### Size Selection

#### **Step 34.**

Air the dry beads for 10 minutes while the tube is on the magnetic stand with the lid open.

#### Size Selection

#### **Step 35.**

Elute the DNA target from the beads into 28  $\mu$ l of 10 mM Tris-HCl or 0.1 X TE, pH 8.0.

#### NOTES

## Isabel Gautreau 26 Jan 2015

Be sure not to transfer any beads. Trace amounts of bead carry over may affect the optimal performance of the polymerase used in the NEBNext High-Fidelity 2X PCR Master Mix in the subsequent PCR step.

#### Size Selection

## **Step 36.**

Mix well on a vortex mixer or by pipetting up and down.

#### Size Selection

## **Step 37.**

Quickly spin the tube and place it on a magnetic stand.

## **O DURATION**

00:05:00

#### NOTES

## Isabel Gautreau 29 Sep 2014

Be sure not to transfer any beads. Trace amounts of bead carry over may affect the optimal performance of the polymerase used in the NEBNext High-Fidelity 2X PCR Master Mix in the subsequent PCR step.

## Size Selection

#### **Step 38.**

After the solution is clear (about 5 minutes), transfer 23 µl to a new PCR tube for amplification.

#### **ANNOTATIONS**

## **Lenny Teytelman** 10 Nov 2015

According to the <u>comments</u>, this is a safe stopping point. After the ligated fragments have been purified, it's ok to freeze overnight and start the PCR amplification step the following day.

#### **PCR** Amplification

## Step 39.

Mix the following components in sterile strip tubes (**Total volume** 50  $\mu$ l):

## **■** AMOUNT

25 μl Additional info:

#### ✓ PROTOCOL

## Mixture for PCR Amplification (E7370)

**CONTACT: New England Biolabs** 

## NOTES

## Isabel Gautreau 26 Jan 2015

The primers are provided in NEBNext Singleplex (<u>NEB #E7350</u>) or Multiplex (<u>NEB #E7335</u>, <u>#E7500</u>) Oligos for Illumina.

## Step 39.1.

Adaptor Ligated DNA Fragments 23 µl

## Step 39.2.

NEBNext High Fidelity 2X PCR Master Mix 25 μl

## **■** AMOUNT

25 μl Additional info:



NEBNext High-Fidelity 2X PCR Master Mix - 50 rxns M0541S by New England Biolabs

## Step 39.3.

Index Primer 1 µl

## Step 39.4.

Universal PCR Primer 1 µl

## **PCR** Amplification

## Step 40.

PCR using the following cycling conditions:

## CYCLE STEP TEMPTIME CYCLES Initial Denaturation98°C 30 seconds1 Denaturation 98°C 10 seconds6-15\*

Denaturation 98°C 10 seconds6-15 Annealing 65°C 30 seconds Extension 72°C 30 seconds Final Extension 72°C 5 minutes 1

Hold 4°C ∞

## NOTES

## Isabel Gautreau 26 Jan 2015

We suggest 6 PCR cycles for 1  $\mu$ g DNA input 10 cycles for 50 ng, and 13–15 for 5 ng DNA input. Further optimization of PCR cycle number may be required.

## Cleanup of PCR Amplification

## Step 41.

Vortex AMPure XP beads to resuspend.

## Cleanup of PCR Amplification

#### Step 42.

Add 50  $\mu$ l of resuspended AMPure XP beads to the PCR reactions ( 50  $\mu$  l).

#### Cleanup of PCR Amplification

## Step 43.

Mix well by pipetting up and down at least 10 times.

## Cleanup of PCR Amplification

#### **Step 44.**

Incubate for 5 minutes at room temperature.

© DURATION 00:05:00

## Cleanup of PCR Amplification

#### Step 45.

Quickly spin the tube and place it on an appropriate magnetic stand to separate beads from supernatant.

© DURATION

00:05:00

## Cleanup of PCR Amplification

## Step 46.

After the solution is clear (about 5 minutes), carefully remove and discard the supernatant. Be careful not to disturb the beads that contain DNA targets (**Caution do not discard beads**).

## Cleanup of PCR Amplification

## **Step 47.**

Wash #1: Add 200 µl of 80% ethanol to the PCR plate while in the magnetic stand.

## Cleanup of PCR Amplification

## **Step 48.**

Wash #1: Incubate at room temperature for 30 seconds.

© DURATION

00:00:30

## Cleanup of PCR Amplification

## Step 49.

Wash #1: Carefully remove and discard the supernatant.

## Cleanup of PCR Amplification

## Step 50.

Wash #2: Add 200 µl of 80% ethanol to the PCR plate while in the magnetic stand.

## Cleanup of PCR Amplification

## Step 51.

Wash #2: Incubate at room temperature for 30 seconds.

**O DURATION** 

00:00:30

## Cleanup of PCR Amplification

#### Step 52.

Wash #2: Carefully remove and discard the supernatant.

## Cleanup of PCR Amplification

## **Step 53.**

Air dry the beads for 10 minutes while the PCR plate is on the magnetic stand with the lid open.

© DURATION

00:10:00

## Cleanup of PCR Amplification

#### Step 54.

Elute DNA target from beads into 33 μl 10 mM Tris-HCl, pH 8.0 or 0.1X TE.

#### Cleanup of PCR Amplification

## Step 55.

Mix well by pipetting up and down at least 10 times.

## Cleanup of PCR Amplification

## **Step 56.**

Quickly spin the tube and place it on an appropriate magnetic stand to separate beads from supernatant.

© DURATION

00:05:00

## Cleanup of PCR Amplification

#### **Step 57.**

After the solution is clear (about 5 minutes), carefully transfer 28 µl supernatant to a new PCR tube.

## Cleanup of PCR Amplification

## Step 58.

Dilute the library 5 fold with water, and check the size distribution on an Agilent high sensitivity chip.

## Warnings

This protocol is suitable for input amounts of 50ng or higher. For lower input amounts we recommend that the protocol without size selection be followed.