

# Primer extension assay and sequencing by dideoxynucleotide incorporation on RNA templates

Matus Valach

## Abstract

The identity of a nucleotide or the presence of a bulky modification or strand break in an RNA can be determined by several approaches. When the 3' region of the analyzed RNA is known, extending a (radio)labeled primer by reverse transcription and analyzing the reaction products using denaturing polyacrylamide gel electrophoresis and subsequent (radio)imaging allows mapping of cDNA chain-termination sites. (Based on the design of the experiment, these can be, for example, a strand break, 5' end, pause-inducing site, or site of dideoxynucleotide incorporation.) This procedure is adapted from the previously published protocol by Walker and Lorsch (DOI: 10.1016/B978-0-12-420037-1.00020-8).

**Citation:** Matus Valach Primer extension assay and sequencing by dideoxynucleotide incorporation on RNA templates. **protocols.io**

dx.doi.org/10.17504/protocols.io.ejtbcnn

**Published:** 23 Feb 2016

## Guidelines

- Estimate the concentration of the target in the sample, whether working with a purified transcript or a complex RNA mixture (e.g. total RNA, poly(A) RNA.). For annealing, RNA substrate and labeled oligo should be equimolar. If primer annealing is problematic because of prominent RNA folding (e.g. if the target is a highly structured RNA), the ratio target:primer can be varied in the range of 1:1 to 1:10 to find the optimal hybridization conditions.

## Before start

- Prepare 10× annealing buffer: 0.5 M HEPES, pH 7.5 + 1 M KCl.
- Prepare 5× extension buffer: 0.5 M Tris, pH 8.0 + 50 mM MgCl<sub>2</sub> + 50 mM DTT.
- Prepare 1 M NaOH and 1M HCl.
- Prepare 2× STOP/loading buffer: 95% deionized formamide, 0.025% SDS, 0.5 mM EDTA, 0.025% bromophenol blue, 0.025% xylene cyanol FF.
- Prepare 5× TBE: 49.5 mM Tris + 49.5 mM boric acid + 2 mM EDTA, pH 8.0.
- Use RNase-free water in all solutions.

## Materials

 Potassium chloride [View](#) by [P212121](#)

- 🦋 T4 Polynucleotide Kinase Reaction Buffer - 4.0 ml [B0201S](#) by [New England Biolabs](#)
- 🦋 AMV Reverse Transcriptase - 200 units [M0277S](#) by [New England Biolabs](#)
- 🦋 T4 Polynucleotide Kinase - 500 units [M0201S](#) by [New England Biolabs](#)
- HEPES [BP310](#) by [Fisher Scientific](#)
- EDTA, disodium salt, dihydrate [S312-500](#) by [Fisher Scientific](#)
- Boric acid BP1681 by [Fisher Scientific](#)
- Sodium Hydroxide BP359500 by [Fisher Scientific](#)
- Tris(hydroxymethyl)aminomethane 327360010 by [Thermo Fisher Scientific](#)
- ✓ Hydrochloric Acid by Contributed by users
- NucAway™ Spin Columns AM10070 by [Thermo Scientific](#)
- Magnesium Chloride AC223210010 by [Fisher Scientific](#)
- Dithiothreitol BP1725 by [Fisher Scientific](#)
- Formamide (deionized) BP228100 by [Fisher Scientific](#)
- Sodium Dodecyl Sulfate BP166 by [Fisher Scientific](#)
- Bromophenol Blue BP11525 by [Fisher Scientific](#)
- Xylene Cyanol FF AC42269-0050 by [Fisher Scientific](#)

## Protocol

### Step 1.

To radioactively label the primer at its 5' end, mix the following components (20 µl):

Component	Amount [µl]	Final concentration
10× polynucleotide kinase (PNK) buffer	2	1×
100 µM DNA oligonucleotide	0.5	2.5 µM
6000 Ci/mmol (10 µCi/µl) [ $\gamma$ - <sup>32</sup> P]ATP	5	2.5 µCi/µl
10 U/µl T4 PNK	1	0.5 U/µl
Water (ddH <sub>2</sub> O)	to 20	

### Step 2.

Incubate for 45 minutes at 37 °C.

### Step 3.

Inactivate the enzyme for 20 minutes at 65 °C.

### Step 4.

Remove unincorporated radioactive ATP with a Sephadex column (e.g. NucAway) according to instructions.

### Step 5.

To anneal the labeled primer to the RNA substrate, mix the following components (7 µl):

Component	Amount	Final concentration
10× annealing buffer	0.7 µl	1×
RNA (target)	0.1-20 pmol	

Labeled DNA oligo	0.1-20 pmol
RNase-free water (ddH <sub>2</sub> O)	to 7 µl

### Step 6.

Denature for 1 minute at 95 °C, then slowly cool the sample (at a rate of 0.1 °C/s) to [primer  $T_m$ -5] °C and incubate for additional 10 minutes (e.g. if the primer  $T_m$  is 55 °C, stop at 50 °C).

### Step 7.

After cooling to 42 °C, immediately add the following mixture (to the final volume of 25 µl):

#### 7.1. For regular primer extension:

Component	Amount [µl]	Final concentration
5× extension buffer	5	1×
10 mM (each) dNTPs	2.5	1 mM (each)
10 U/µl AMV reverse transcriptase	1	0.4 U/µl
RNase-free water (ddH <sub>2</sub> O)	9.5	

#### 7.2. For dideoxy-termination sequencing (e.g. A lane):

Component	Amount [µl]	Final concentration
5× extension buffer	5	1×
10 mM (each) dNTPs	2.5	1 mM (each)
5 mM ddTTP	1.25	0.25 mM
10 U/µl AMV reverse transcriptase	1	0.4 U/µl
RNase-free water (ddH <sub>2</sub> O)	8.25	

### Step 8.

Incubate for 45 minutes at 42 °C.

### Step 9.

Add 3 µl of 1 M NaOH (final concentration 100 mM) to degrade RNA.

### Step 10.

Incubate for 10 minutes at 95 °C.

### Step 11.

Add 3 µl of 1 M HCl to neutralize the solution.

### Step 12.

Add 31 µl of 2× STOP/loading buffer.

### Step 13.

Incubate for 5-10 minutes at 95 °C.

### Step 14.

Place on ice for 3 minutes.

### Step 15.

Load on a pre-ran denaturing PAGE gel (1× TBE, 7M urea), run the electrophoresis, and then expose the gel to a film or phosphor-imaging screen.

## Warnings

When working with radioactive material, follow the proper radiation safety procedures.