

# PCR on Euplotes whole cells

Angela Piersanti

## Abstract

**Citation:** Angela Piersanti PCR on Euplotes whole cells. **protocols.io**

dx.doi.org/10.17504/protocols.io.hhvb366

**Published:** 01 Apr 2017

## Protocol

### Step 1.

Collect ~ 2 ml of cells from a culture at 400 cells/ml, concentrate them by centrifugation, 7 min at 3000 rpm. Remove the supernatant by pipetting.

### Step 2.

Wash cells twice by adding ~ 1 ml of ddH<sub>2</sub>O. Centrifuge for 7 min at 3000 rpm, remove the supernatant by pipetting.

### Step 3.

Put 15 µl (~20 *Euplotes* cells) in a PCR tube and put at 98°C for 20 min in the thermal cycler.

### Step 4.

After that add the PCR mix to the cells, 50 µl of final volume for *Euplotes crassus* and 20 µl of final volume for *Euplotes focardii*.

### Step 5.

Add the polymerase and proceed with the thermal cycle.