

## Agarose gel electrophoresis

Cristian Riccio<sup>1</sup>

<sup>1</sup>University of Cambridge

[dx.doi.org/10.17504/protocols.io.3ddgi26](https://doi.org/10.17504/protocols.io.3ddgi26)

Working



Cristian Riccio ⚡

### ABSTRACT

Agarose gel for DNA electrophoresis

This protocol is useful to check the approximate size of a mixture of DNA fragments.

For example, after a HMW gDNA extraction, it is good to check that the DNA we extracted is actually high-molecular weight. Another use is for separating PCR amplicons when checking for a mutation. Primers can be designed in such a way that PCR will yield two fragments of different sizes depending on whether a sequence is present or not. If a deletion is present, then the amplified fragment will be smaller than for the DNA where there is no deletion.


### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Gel Loading Dye Orange (6X) - 4.0 ml	B7022S	New England Biolabs
TAE Buffer (Tris-acetate-EDTA)	B49	
Hyperladder 4 (100bp)	BIO-33056	Bioline
Agarose		Sigma
Image Lab Software	View	Bio-rad Laboratories
Autoclave Indicator Tapes	AT011.SIZE.1ROLLS	Bio Basic Inc.
Ethidium bromide 10 mg/ml	E1510	Sigma Aldrich
PowerPac™ HC Power Supply	1645052	BioRad Sciences
Hyperladder 1kb	BIO33025	Bioline
agarose gel container	View	
Bio-Rad electrophoresis machine	View	
Bio-Rad Gel Doc XR with Image Lab Software	View	
Microwave oven	View	

### SAFETY WARNINGS

Ethidium bromide is a carcinogen and you should be trained on how to handle it before use. EtBr also needs to be disposed of in a special way.

- 1 Get a 250 ml flask and a 250 ml measuring cylinder.

- 2 Add 1 gram agarose to the flask. Add 130 ml TAE in measuring cylinder. Taper container. Put comb in container.
- 3 Microwave flask with agarose + TAE for 1 minute. After one minute, check that the solution is homogeneous. If not, microwave in steps of 15 seconds until it is homogeneous.
- 4 Move taped and combed container to the cold room along with the quasi-boiling agarose solution. Add 1 microlitre ethidium bromide to the agarose solution. Pour the agarose solution into the container. Make sure there is no spillage. If there is spillage, stop pouring, clean the mess and start again.
- 5 Add a small amount of water to the flask, rinse it and dispose of the contaminated water in the liquid EtBr waste through the white material.
- 6 Wash flask with water and leave it in the dirty dishes trolley.
- 7 Leave the gel to solidify in the cold room. If you want to use the gel the following day, wrap it in a plastic sheet to protect it from desiccation. 2h
- 8 Remove tape, put gel in electrophoresis tank. Add 5 microlitres ladder.
- 9 Add 2 microlitres loading buffer to 5 microlitres of DNA sample (~ 100 nanograms) and load into gel.
- 10 Attach the cables from the tank to the Bio-Rad power machine. Set the programme to 150 V, 400 mA and 30 minutes. Start run. 30m
- 11 Check gel every 10 minutes. Put gel in Imager. Open Image Lab. New Protocol, Position gel, zoom in/out. Start protocol. Image is ready on computer. Invert contrast, change contrast, print image. Don't touch computer without gloves. It is contaminated.  
 **Bio-Rad 2019-05-26 15hr 05min.pdf**

#### Acknowledgments

- 12 Jérémie Le Pen and other colleagues in the Miska lab  
Ulrika Yuan for recording



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited