

# Formulation of MBBM (Modified Bold's Basal Medium)

David Dunigan and Irina Agarkova

## Abstract

This medium is for growing *Chlorella variabilis* NC64A, *Chlorella heliozoae* SAG 3.83, and *Coccomyxa* C-169.

**Citation:** David Dunigan and Irina Agarkova Formulation of MBBM (Modified Bold's Basal Medium). **protocols.io**  
dx.doi.org/10.17504/protocols.io.etwbepe

**Published:** 13 Jun 2016

## Guidelines

### BBM Stock Solutions

1. 25.0 gm  $\text{NaNO}_3$  per liter d- $\text{H}_2\text{O}$
2. 2.5 gm  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$  per liter d- $\text{H}_2\text{O}$
3. 7.5 gm  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  per liter d- $\text{H}_2\text{O}$
4. 7.5 gm  $\text{K}_2\text{HPO}_4$  per liter d- $\text{H}_2\text{O}$
5. 17.5 gm  $\text{KH}_2\text{PO}_4$  per liter d- $\text{H}_2\text{O}$
6. 2.5 gm  $\text{NaCl}$  per liter d- $\text{H}_2\text{O}$
7. 50.0 gm disodium EDTA, 31.0 gm  $\text{KOH}$  per liter d- $\text{H}_2\text{O}$
8. 4.98 gm  $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$  per liter acidified  $\text{H}_2\text{O}$  (Acidified  $\text{H}_2\text{O}$  is 999.0 mL d- $\text{H}_2\text{O}$  + 1.0 mL concentrated  $\text{H}_2\text{SO}_4$ )
9. 11.42 gm  $\text{H}_3\text{BO}_3$  per liter d- $\text{H}_2\text{O}$
10. 8.82 gm  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ , 1.44 gm  $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ , 0.71 gm  $\text{MoO}_3$ , 1.57 gm  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ , and 0.49 gm  $\text{Co}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$  per liter d- $\text{H}_2\text{O}$

Note: Stock solution 10 takes weeks for all of the salts to dissolve. Use as suspension until then.

### BBM Preparation

to 950 mL of d- $\text{H}_2\text{O}$  add:

10.0 mL of stock solutions 1, 2, 3, 4, 5 and 6

1.0 mL of stock solutions 7, 8 and 9

2.0 mL of stock solution 10

### MBBM Preparation

to 950 mL of d- $\text{H}_2\text{O}$  add:

10.0 mL of stock solutions 1, 2, 3, 4, 5 and 6

1.0 mL of stock solutions 7, 8 and 9

2.0 mL of stock solution 10

1.0 gm of bacto-peptone

5.0 gm of sucrose

Tetracycline (filter sterilized, 10 µg/mL final concentration) is added after the media is autoclaved and cool.

For MBBM plates, agar is added to 1.5% before autoclaving.

For MBBM soft agar (for titering), agar is added to 0.75% before autoclaving.

## Protocol