



Oct 14, 2019

## Western Blot (Tank-Blot)

[iGEM Dusseldorf<sup>1</sup>](#)<sup>1</sup>Heinrich-Heine Universität Düsseldorf

1

Works for me

[dx.doi.org/10.17504/protocols.io.765hr6](https://dx.doi.org/10.17504/protocols.io.765hr6)

iGEM Dusseldorf ⚡

### ABSTRACT

#### Blotting Buffer:

- 100mL 10x SDS running buffer
  - 200mL 100% Methanol
  - 700mL H<sub>2</sub>O
- 
- Fill a tray with blotting buffer
  - Transfer SDS gel to the blotting buffer
  - Assemble blotting sandwich (take it, soak sponge and place on it, soak filter paper and place on it, place inverted SDS gel on top, without creating bubbles. Add charged PDVF membrane with forceps ( charged in 100% ethanol and then soaked in buffer )
  - Add soaked filter
  - Streak out bubbles
  - Add soaked sponge
  - Press together and close container
  - Place in blotting apparatus and add blotting buffer
  - Run for 1h at 100V 350-500mA



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited