

# ষ্ট Earth Microbiome Project (EMP) DNA extraction protocol

Lisa Marotz, Tara Schwartz, Luke Thompson, Greg Humphrey, Grant Gogul, James Gaffney, Amnon Amir, Rob Knight

#### **Abstract**

Citation: Lisa Marotz, Tara Schwartz, Luke Thompson, Greg Humphrey, Grant Gogul, James Gaffney, Amnon Amir, Rob

Knight Earth Microbiome Project (EMP) DNA extraction protocol. protocols.io

dx.doi.org/10.17504/protocols.io.pcydixw

Published: 06 Aug 2018

# **Guidelines**

# MoBio PowerMag Soil DNA Isolation Kit (Optimized for KingFisher)

The Knight lab has transitioned to the PowerMag Soil DNA Isolation Kit (Optimized for KingFisher). We have validated a variety of sample types to ensure reproducibility when compared to MoBio PowerSoil Extraction Kit. This transition occurred to increase efficiency and reduce DNA extraction time from 6-8 hours to 2.5-3 hours.

The protocol is followed as MoBio recommends, with an added 10-minute water bath at 65°C after step 7.

The new kit can be implemented on the epMotion using a magnetic bead plate adapter for the epMotion. However, it does not reduce the amount of time the extraction takes by a significant amount: tests showed that the new kit on the epMotion took longer than the old kit. Comparison showed both the old and new kit performed well on the epMotion.

#### PowerSoil-htp 96 Well Soil DNA Isolation Kit

The MoBio PowerSoil DNA Isolation Kit is still compatible if KingFisher instrumentation is not available.

# Items included in the extraction kit (for 1-96 well extraction)

- (1) Bead plate
- (1) Spin plate (filter)
- (1) 0.5ml collection
- (4) 1.0ml collection plates
- (2) 2.0ml collection plates
- (1) Microplate (DNA elution)
- Sealing Tape
- Centrifuge Tape
- Elution Sealing Mat
- Labeled solutions

# **Important Considerations**

- Normal diameter 1 ml pipet tips are too large for some of the pipetting steps. To get around this
  problem we use a Rainin 300 ul 8-channel with filtered tips (Rainin #SR-L300F) and a Rainin
  1000 ul 8-channel with extended length filtered tips (Rainin #RTS-L1000XF).
- 2. Make sure that all of the necessary consumables and reagents are in place before you start the extraction. Remember, each pipetting step will require 1 box of 96 tips per plate.
- 3. We use individually wrapped reagent reservoirs and expose them to UV light for 30 minutes prior to usage.
- 4. In Step 10, it is important that the plate not rub against any surfaces in the shaker.
- 5. Make sure that the alpha-numeric grid is in the same orientation across all the plates. On a few occasions we have noticed that the sticker is not always in the same position on the plate.

#### **Next Steps**

- 16S rRNA Amplification Protocol
- 18S rRNA Amplification Protocol

#### **Before start**

- 1. Clean all surfaces and pipettors to remove DNA
- 2. Label Plates
  - ∘ Plate #1 1ml collection plate
  - ∘ Plate #2 1ml collection plate
  - Plate #3 1ml collection plate
  - ∘ Plate #4 1ml collection plate
  - Plate #5 2ml collection plate
  - ∘ Plate #6 2ml collection plate
- 3. UV Sterilize and Label reservoirs
  - Bead solution 750ul/well
  - o C1 60ul/well
  - o C2 250ul/well

- o C3 200ul/well
- C4 650ul/well
- C5 500ul/well
- o C6 100ul/well

#### **Protocol**

# Step 1.

BEFORE THE FIRST USE ONLY, Solution C5-D must be prepared. Add an equal amount of 100% Ethanol to Solution C5-D (for the 4 prep kit = 120 ml, or for the 12 prep kit = 360 ml). Mix well. Put a check mark in the "ethanol added" box on the bottle cap label.

# Step 2.

Centrifuge Bead Plate for 1 min at 2500 x g to pellet the beads.

**O DURATION** 

00:01:00 : Centrifugation

Step 3.

Remove the Square Well Mat from the PowerSoil®-htp Bead Plate and set aside.

#### Step 4.

Add 0.1 to 0.25 grams of soil sample or sample swab.

#### NOTES

This is an appropriate stopping point and you can store the PowerSoil®-htp Bead Plate at 4°C covered with the Square Well Mat. This is the most time consuming step of the protocol. Care must be taken to avoid cross contamination between sample wells.

# Step 5.

Add 750 µl of PowerSoil®-htp Bead Solution to the wells of the PowerSoil®-htp Bead Plate.

**■** AMOUNT

750 µl: PowerSoil®-htp Bead Solution

#### Step 6.

Check Solution C1. If Solution C1 has precipitated, heat solution at 60°C until the precipitate has dissolved.

**↓** TEMPERATURE

60 °C: Heating solution

NOTES

Solution C1 contains SDS. If it gets cold, it will precipitate. Heating at 60°C will dissolve the SDS. Solution C1 can be used while it is still warm.

### Step 7.

Add 60 µl of Solution C1. Secure the Square Well Mat (from step 3) tightly to the plate.

**■** AMOUNT

60 µl: Solution C1

# Step 8.

Place sealed plates in 65°C water bath for 10 min. DO NOT SUBMERGE THE PLATES.

© DURATION

00:10:00 : Water bath

# Step 9.

Place PowerSoil®-htp Bead Plate between the aluminum plate adapters and securely fasten to the 96 Well Plate Shaker.

#### Step 10.

Shake at speed 20 for 20 minutes.

© DURATION

00:20:00 : Shaking

NOTES

It is important that the plate not rub against any surfaces in the shaker.

#### Step 11.

Centrifuge at room temperature for 6 minutes at 4500 x g. While centrifuging, aliquot 250  $\mu$ l of Solution C2 into each well of Plate #1 and cover with Sealing Tape.

**■** AMOUNT

250 μl : Solution C2

© DURATION

00:06:00 : Centrifugation

#### NOTES

The Sealing Tape can be re-used when centrifuging Plate #1 in step 14 if handled carefully.

# Step 12.

Remove and discard the Square Well Mat from the Bead Plate.

# Step 13.

Carefully remove the Sealing Tape from Plate #1 and transfer the supernatant ( 400-500µl) from the Bead Plate to Plate #1 and pipette up and down 4 times.

#### NOTES

The supernatant may still contain some particles.

# Step 14.

Re-apply the Sealing Tape to Plate #1. Incubate at 4°C for 10 minutes.

↓ TEMPERATURE4 °C : IncubationĠ DURATION

00:10:00 : Incubation

# Step 15.

Centrifuge Plate #1 at room temperature for 6 minutes at 4500 x g. While centrifuging, aliquot 200  $\mu$ l Solution C3 into each well of Plate #3, then cover with Sealing Tape.

# **■** AMOUNT

200  $\mu l$  : Solution C3

© DURATION

00:06:00 : Centrifugation

# **P** NOTES

The Sealing Tape can be re-used when centrifuging Plate #3 in step 16 if handled carefully.

#### Step 16.

After centrifugation, carefully remove and discard Sealing Tape from Plate #1.

#### Step 17.

Avoiding the pellet, transfer the entire volume (600  $\mu$ l depending on sample type) of supernatant in Plate #1 to Plate #2.

# Step 18.

Apply Sealing Tape to Plate #2 and centrifuge at room temperature for 6 minutes at 4500 x g.

#### **O** DURATION

00:06:00: Centrifugation

#### Step 19.

Carefully remove Sealing Tape from Plate #2 and Plate #3.

#### Step 20.

Avoiding the pellet, transfer the entire volume of supernatant (600  $\mu$ l) from Plate #2 to Plate #3 and pipette up and down 4 times.

# Step 21.

Re-apply Sealing Tape to Plate #3. Incubate at 4°C for 10 minutes.

↓ TEMPERATURE4 °C : IncubationĠ DURATION

00:10:00 : Incubation

Step 22.

Centrifuge at room temperature for 6 minutes at 4500 x g.

**O DURATION** 

00:06:00 : Centrifugation

Step 23.

Carefully remove and discard Sealing Tape from Plate #3.

#### Step 24.

Avoiding the pellet, transfer the entire volume of supernatant (750 µl) to Plate #4.

# Step 25.

Apply Sealing Tape to Plate #4 and centrifuge at room temperature for 6 minutes at 4500 x g. While centrifuging, add 650  $\mu$ l of Solution C4 to Plate #5.

**■** AMOUNT

650 μl: Solution C4

**O** DURATION

00:06:00 : Centrifugation

Step 26.

Avoiding any residual pellet, transfer up to 650 µl of supernatant in Plate #4 to Plate #5.

#### Step 27.

Add a second 650 µl (1300 µl C4 total) aliquot of Solution C4 to each well of Plate #5.

**■** AMOUNT

650 µl: Solution C4

NOTES

It is safe to stop the protocol at this step and store the samples covered with Sealing Tape at 4°C.

Make sure to briefly centrifuge the plate to collect any condensate on the plate seal after overnight storage.

### Step 28.

Pipet samples "up and down" to mix.

#### Step 29.

Place Spin Plate onto Plate #6.

#### Step 30.

Load approximately 650 µl from Plate #5 into each well of the Spin Plate and apply Centrifuge Tape.

#### Step 31.

Centrifuge at room temperature for 5 minutes at 4500 x g.

© DURATION

00:05:00: Centrifugation

#### Step 32.

Discard the flow through and place the Spin Plate back on Plate #6. Carefully remove and discard the Centrifuge Tape.

#### Step 33.

Repeat steps 30-32 until all the supernatant has been processed. Discard the final flow through.



Repeating steps 30-32 -> go to step #30

#### Step 34.

Place the Spin Plate back on Plate #6.

# Step 35.

Confirm that ethanol has been added to Solution C5-D (see step 1). Add 500  $\mu$ l of Solution C5-D to each well of the Spin Plate. Apply Centrifuge Tape to the Spin Plate.



500 μl: Solution C5-D

# Step 36.

Centrifuge at room temperature for 5 minutes at 4500 x g.

**O DURATION** 

00:05:00 : Centrifugation

#### Step 37.

Discard the flow through and place the Spin Plate back on Plate #6.

#### Step 38.

Centrifuge again at room temperature for 6 minutes at 4500 x g.

#### **O** DURATION

00:06:00 : Centrifugation

# Step 39.

Discard the flow through.

# Step 40.

Carefully place the Spin Plate onto the Microplate. Remove Centrifuge Tape from the Spin Plate and discard.

#### Step 41.

Add 100 µl of Solution C6 to the center of each well of the Spin Plate. Apply Centrifuge Tape.

# **■** AMOUNT

100 µl: Solution C6

#### Step 42.

Let C6 sit on the filter for 10 minutes at room temperature before final centrifugation step.

#### **O** DURATION

00:10:00 : C6 sit on filter

#### Step 43.

Centrifuge at room temperature for 7 minutes at 4500 x g.

#### **O** DURATION

00:07:00 : Centrifugation

# Step 44.

Remove Centrifuge Tape and discard.

#### Step 45.

Cover wells of Microplate with the Elution Sealing Mat provided. DNA is now ready for any downstream application. No further steps are required.

#### NOTES

Prolonged storage at 4°C will result in the evaporation of eluted DNA. We recommend storing DNA

frozen (-20°C or -80°C). Solution C6 does not contain EDTA. To concentrate the DNA see the Hints and Troubleshooting Guide provided in the MoBio protocol.

# **Warnings**

Please wear gloves at all times.

Please refer to the SDS (Safety Data Sheet) for hazard information.