



Sep 04, 2019

Top2 inhibitor sensitivity

Jacob Kirkland¹

¹Stanford University

1

Works for me

dx.doi.org/10.17504/protocols.io.k5zcy76



Jacob G. Kirkland
Stanford University



ABSTRACT

Protocol for Top2i sensitivity in attached cells in 96-well format. Short term exposure to Top2i (~3 days)

MATERIALS

NAME	CATALOG #	VENDOR
50 mg Doxorubicin HCl	orb340472	biorbyt
500 mg Etoposide	orb322762	biorbyt
Doxycycline Hyclate	D-500	Gold Biotechnology

SAFETY WARNINGS

TRPz lenti-virus is BSL2

BEFORE STARTING

Using a Tripz doxyrubicin inducible shRNA system that is already in transduced cells.

- 1 Plate (2) 6cm dishes at 25% confluency (approx. 500K cells). Let attach overnight

Induction of shRNAs

- 2 Add 1000x Doxycyclin to 1 plate of cells with fresh media and induce for 72 hours changing media with fresh doxycyclin daily. Also change media to corresponding non-induced plate

🕒 72:00:00

In 96-well plates

- 3 Split cells at density appropriate for cell type. (eg. 5000k cells/well) in total of 150ul of media per well. Make master mixes of cells/well in 50ml conicals and distribute using multi-channel pipette. Add doxycyclin to induced samples at 500x. Let cells attach overnight

🕒 14:00:00

14h

- 4
 - Dilute Top2i in media 96-well round bottom plate for transfer to plates with cells
 - Dilution in media without Doxycycline so that both cell samples can get Top2i from the same dilutions
 - In our case 2ul of stock concentration of drug (see making stock plate) in 250ul of media (1:125) --> (1:500 when 50uls is added to 150ul already in well.
 - Doxycycline concentration will be diluted in induced samples but was starting with a higher concentration).
 - One row is enough for 1/2 plate. Make multiple rows of dilutions for more plates.

5 ⌚ 48:00:00

Incubate cells with Top2i for 2 days

- 6
- Dilute cell titer blue 1:2 in PBS
 - Add 20ul of cell titer blue to wells and let sit overnight

⌚ 14:00:00

- 7
- Transfer supernatant to 384-well black plates for reading
 - Excite at 544nm and emission of 590nm



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited