

Clonal virus purification

Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen

Abstract

For use in "Obtaining pure cyanophage stocks (liquid assay)"

Citation: Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen Clonal virus purification. **protocols.io**
dx.doi.org/10.17504/protocols.io.dqp5vm

Published: 08 Feb 2016

Protocol

Step 1.

Use 13-x-100-mm culture tubes (or 24-well plates).

Step 2.

Prepare exponentially growing target cells (e.g., >100 mL).

Step 3.

Dilute some of the titered lysate to 1 infective virus/mL.

Step 4.

Add 0.2 mL (0.2 infectious units) to each of 20 tubes of susceptible host cells.

Step 5.

Monitor tubes for 2 to 3 weeks.

Step 6.

Cultures in which lysis occurs are assumed to be the result of a single-virus infection.

📌 NOTES

Amy Chan 02 Sep 2015

The probability that more than 1 infective unit occurred in a given culture is 0.0176.

Step 7.

If lysis occurs in 4 tubes or less of 20, it is assumed that lysis in each tube was caused by one infectious unit, therefore each tube would contain a separate phage clone.

Step 8.

Propagate an aliquot from all the tubes to confirm the results.

Step 9.

If lysis occurs in more than 4 tubes, repeat the clone out procedure by reducing the volume of diluted lysate added to the 20 tubes.

📌 NOTES

Amy Chan 05 Oct 2015

Add 0.1 mL instead of 0.2 mL.

Step 10.

Scale up each phage clone to make primary phage stocks.

📌 NOTES

Amy Chan 05 Oct 2015

e.g., add 5 μ L of the lysate to 40 mL of cells.

Step 11.

Centrifuge, filter, and titer the stock, store at 4°C in the dark.