

Stacking gel buffer (4X)

DAVID DUNIGAN AND IRINA AGARKOVA

Abstract

For use in [Polyacrylamide Gel System For Electrophoresis Of Proteins](#).

Citation: DAVID DUNIGAN AND IRINA AGARKOVA Stacking gel buffer (4X). **protocols.io**

dx.doi.org/10.17504/protocols.io.eq4bdyw

Published: 13 Jun 2016

Protocol

Step 1.

Dissolve 6.1 gm of Tris base in 90 mL of d-H₂O.

Step 2.

Adjust the pH to 6.8 with concentrated HCl.

Step 3.

Adjust to a final volume of 100 mL.