



Adult mouse skin dissociation protocol (on ice)

Andrew Potter¹

¹CCHMC

Sep 21, 2018 dx.doi.org/10.17504/protocols.io.tupenvn

Working

Human Cell Atlas Method Development Community



Andrew Potter   

ABSTRACT

This protocol was developed to dissociate adult (8-10 wk) mouse skin "on ice". It utilizes two layers of digestion with a *Bacillus licheniformis* protease cocktail, combined with mechanical disruption from a dounce homogenizer. The cell yield is 3000 cells/mg.



Skin adult mouse
dissociation
protocol.pdf

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

***Bacillus Licheniformis* Enzyme Mix (1 mL per 23 mg tissue):**

100 μ L b. lich 100 mg/mL (10 mg/mL final) (Sigma, P5380)

1 μ L 0.5 M EDTA (Sigma, A8806)

899 μ L DPBS (no Ca, Mg) ThermoFisher (cat. #14190)

Preparing enzymes:

The enzyme is made up in DPBS (#14190). It is aliquoted and stored at -80 °C at 100 mg/mL in 100 μ L aliquots..

Reagents

Enzymes, trypsin inhibitor, BSA and DNase are made up in DPBS (no Ca, no Mg) from Thermo Fisher (14190).

Bovine Serum Albumin - Sigma (A8806).

Hypothermosol FRS

Required supplies:

2 mL dounce homogenizer – Bellco (1984-10002)

Centrifuge for 1.5 mL, 15 mL conicals

Pipettes and pipet tips

15 ml Conicals (MLS)

1.5 mL tubes (MLS)

30 μ M filters - Miltenyi (130-098-458)

Petri dishes (MLS)

Razor blades (MLS)

Ice bucket w/ice

Hemocytometers - InCyto Neubauer Improved (DHC-NO1-5)

Isolating tissue

1. After euthanizing mouse, remove hair using Nair: dab with Nair, wait 30 secs, wipe with wet paper towel.
2. Isolate tissue and place in ice-cold hypothermosol.
3. Scrape off underlying layer of fatty / connective tissue using scalpel before proceeding.
4. Mince skin tissue thoroughly on petri dish on ice for 3-4 min on ice into 1-mm³ pieces using razor blade while manipulating tissue with forceps – you will need to use grinding motion and vigorously break up tissue.

 00:04:00 mincing


1st digest layer

5. Place 23 mg minced tissue into 1 mL B. Lich enzyme cocktail. Incubate on ice.

 23 mg minced tissue

6. Shake every min; triturate 10x every 2 min with p1000 w/tip cut (start triturating at 2 min) for 20 min.

 00:20:00 digest

 00:01:00 shake

 00:02:00 triturate 10X

7. After 20 mins of triturating on ice, use pipet to transfer digest mix to 2 mL dounce homogenizer. Use 10 strokes of Pestle A every 2 min (4 series total, 8 min). Digest mix should become turbid.

 00:02:00 dounce homogenize

 00:08:00 digest using dounce

8. Transfer back to 1.5 mL tube using 1 mL serological pipet. Mix thoroughly and allow to settle on ice 2 min.

 00:02:00 settle on ice

9. Save 70% (700 µL) of supernatant, leaving chunks at the bottom of the tube; apply to 30 µM filter on 15 mL conical. Rinse filter w/5 mL ice-cold PBS/BSA 0.04%. Save flow through on ice and keep filter on tube for 2nd layer.

 700 µl save supernatant

 5 ml rinse filter w/ice-cold PBS/BSA 0.04%

2nd digest layer

10. Add additional 1 mL b. Lich enzyme mix to residual tissue chunks.

 1 ml b. lich enzyme mix

11. Triturate 10x every 2 min, shake every min while incubating on ice for 20 additional mins. (50 min. total digest time).

11

🕒 00:20:00 additional digest time

🕒 00:02:00 triturate 10x

🕒 00:01:00 shake every min.

12

12. Transfer entire volume to same 30 μ M filter on 15 mL conical. Rinse with additional 5 mL ice-cold PBS/BSA 0.04%.

🧊 5 ml ice-cold PBS/BSA 0.04%

Preparing cells for single cell analysis

13

13. Centrifuge at 300 g for 5 min at 4 °C. Remove supernatant & re-suspend in 100 μ L PBS/BSA 0.04%. Examine using hemocytometer with trypan blue.

🕒 00:05:00 centrifuge at 300 g

🧊 100 μ L re-suspend in ice-cold PBS/BSA 0.04%

14

14. Adjust concentration to 1,000 cells/ μ L for Chromium or 100 cells/ μ L for DropSeq.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited