

M9 Minimal Media

Dr. Steven Wilhelm

Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol

Modified from Pardee, A. B., F. Jacob, and J. Monod. 1959. The genetic control and cytoplasmic expression of "inducibility" in the synthesis of β -galactosidase in *E. coli*. *J. Mol. Biol.* 1:165-178

Citation: Dr. Steven Wilhelm M9 Minimal Media. **protocols.io**

dx.doi.org/10.17504/protocols.io.ibqcamw

Published: 06 Jun 2017

Protocol

Step 1.

Add 200 mL 5x M9 salts to a clean media bottle

PROTOCOL

. [M9 Salts](#)

CONTACT: [Steven Wilhelm](#)

Step 1.1.

Add 1 L dH₂O to a clean media bottle

Step 1.2.

Add 64 g Na₂HPO₄·7H₂O


REAGENTS

 Disodium hydrogen phosphate by Contributed by users

Step 1.3.

Add 15 g KH₂PO₄

REAGENTS

 Potassium phosphate (monobasic) [View](#) by [P212121](#)

Step 1.4.

Add 2.5 g NaCl



REAGENTS

✓ Sodium Chloride [PubChem CID: 5234](#) by Contributed by users

Step 1.5.

Add 5 g NH₄Cl



REAGENTS

✓ Ammonium Chloride [View](#) by [P212121](#)

Step 1.6.

Autoclave at 121°C for 20 m

Step 2.

Add 2 mL 1 M MgSO₄



REAGENTS

✓ Magnesium sulfate heptahydrate by Contributed by users

Step 3.

Add 1 mL 0.1 M CaCl₂



REAGENTS

✓ Calcium Chloride by Contributed by users

Step 4.

Add sterile dH₂O to bring solution to 1 L final volume

Step 5.

Autoclave at 121°C for 20 m