

CZV8- fungal protoplast agar

Oliver Mead

Abstract

CZV8-proto agar

	500 mL	300 mL
Czapek-Dox broth	22.7 g	13.62 g
agar	5 g	3 g
sorbitol	91.1 g	54.7 g
centrifuged V8 juice	100 mL	60 mL
add H ₂ O to	500 mL	300 mL
adjust to pH 6.0 with NaOH		
autoclave and pour thin plates (approx. 10mL)		


Citation: Oliver Mead CZV8- fungal protoplast agar. **protocols.io**

dx.doi.org/10.17504/protocols.io.f7dbri6

Published: 25 Oct 2016

Materials

Czapek Dox [CM0097](#) by [Thermo Scientific](#)

 Campbells V8 vegetable Juice (original) by Contributed by users

Protocol

Step 1.