

Binding Buffer (3M GuHCl, 3.75M NH₄Ac pH 6)

Sarah Hessen-Schmidt

Abstract

Citation: Sarah Hessen-Schmidt Binding Buffer (3M GuHCl, 3.75M NH₄Ac pH 6). **protocols.io**

dx.doi.org/10.17504/protocols.io.fjsbkne

Published: 10 Aug 2016

Protocol

Step 1.

Add 57.318g Guanidine Hydrochloride to an autoclaved bottle



REAGENTS

✓ Guanidine Hydrochloride GU103 by Contributed by users

Step 2.

Add 50ml MilliQ

Step 3.

Add 50ml 7.5M Ammonium Acetate



PROTOCOL

. [7.5M Ammonium Acetate](#)

CONTACT: [Sarah Hessen-Schmidt](#)

Step 3.1.

57.81g Ammonium Acetate (NH₄OAc) in 100ml H₂O



AMOUNT

58 g Additional info:



REAGENTS

Ammonium Acetate A1542-500G by [Sigma Aldrich](#)