

Whole protein extraction from tissues

Vivian Liu

Abstract

This protocol is for whole protein extraction from tissues like kidney, liver, muscle etc. Histone proteins which are tightly twisted on genomic DNA can also be extracted with this protocol.

Citation: Vivian Liu Whole protein extraction from tissues. **protocols.io**

dx.doi.org/10.17504/protocols.io.dkx4xm

Published: 06 Aug 2015

Before start

ULTRA-TURRAX Dispenser T10 basis S25 (IKA)

Rotating mixer

Bench top centrifuge

Protocol

Step 1.

Pre-cool lysis buffer on ice. Add protease inhibitors just before using.

Step 2.

Add 10 fold lysis buffer to tissues (1ml lysis buffer / 100mg tissue)

NOTES

Vivian Liu 06 Aug 2015

[RIPA lysis buffer \(strong\) recipe](#)

Step 3.

Homogenize tissues with dispenser (working time 10s / rest time 20s, repeat 5 times) until no visible tissue blocks.

Step 4.

Mix samples on a rotating mixer at 4 °C for 2h.

Step 5.

Spin @ 16,000g, 10min, 4°C. If necessary, centrifuge again.

Step 6.

Save supernatant. Determine protein conc. .

Step 7.

Make aliquots and store in -80°C.